



650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

Classic Line, Four Drawer, Emergency Cart, Specialty Package #6301



**Color Shown is Red.
Additional Colors Available**

Assembled Cart Dimensions:

- 39.5" H x 35" W x 22" D
(100.3 cm H x 88.9 cm W x 55.8 cm D)
- Cart net weight: 128 lbs (58.0 kg)
- Approximate shipping weight: 167 lbs (75.8 kg)
- Approximate shipping dims: 47" H x 33" W x 26" D
(119.4 cm H x 83.8 cm W x 66 cm D)

Specifications

Construction:

- Sturdy, uni-body steel construction

Features:

- Painted steel break-away lock bar
- Pack of 100 break-away seals
- 5" (12.7 cm) ball-bearing premium casters, all full swivel-one tracking, one braking
- 24" of vertical drawer space
- Cabinet and drawer design allows for customized drawer configurations
- Painted steel top rail / push handle
- Removable plastic top
- Ball-bearing full extension drawer slides
- Stabilizing pontoon bumpers
- Pre-drilled threaded holes for side mounted accessories
- Durable hard-baked powder coat finish
- Cart is delivered mostly assembled
- Twelve year warranty on manufactured parts
- Latex free

Accessory Package Includes:

- Cardiac board and brackets (680406)
- Drop Shelf
- Utility hooks (4 ea) (680402)
- IV pole (680401)
- Oxygen tank holder - D/E size cylinders (AL680403)
- 6 compartment drawer trays (3 ea)
(fit three per drawer), (40108)
- 7 compartment drawer trays (3 ea)
(fit three per drawer), (40109)

Drawer Configuration:

- (2) – 3" (7.6 cm)
- (1) – 6" (15.2 cm)
- (1) – 12" (30.5 cm)

Internal Drawer Dimensions:

- 23" W x 17" D (58.4 cm W x 43.2 cm D)