



650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

Classic Line, Five Drawer, Emergency Cart, with bottom compartment, Specialty Package #6403



Color Shown is Red.
Additional Colors Available

Assembled Cart Dimensions:

- 45.75" H x 35" W x 22" D
(116.2 cm H x 89 cm W x 56 cm D)
- With IV pole: 86.25" (219 cm H)
- Cart weight: 150 lbs (68 kg)
- Approximate shipping weight: 189 lbs (85.7 kg)
- Approximate shipping dims: 53" H x 33" W x 26" D
(134.6 cm H x 84 cm W x 66 cm D)

Specifications

Construction:

- Sturdy, uni-body steel construction

Features:

- Painted steel break-away lock bar
- Pack of 100 break-away seals
- 30" of vertical drawer space
- Cabinet and drawer design allows for customized drawer configurations
- 12" flip-up door on bottom compartment
- Painted steel top rail
- Stabilizing pontoon bumpers
- 5" (12.7 cm) ball-bearing premium casters, all full swivel - one tracking, one with brake
- Removable plastic top
- Ball-bearing full extension drawer slides
- Cart is delivered mostly assembled, IV pole & O2 holder are customer installed
- Pre-drilled threaded holes for side mounted accessories
- Durable, hard-baked powder coat finish
- Latex free
- Twelve year warranty on manufactured parts

Accessory Package Includes:

- Cardiac board and brackets (680406)
- Drop shelf (680302)
- Utility hooks (4 ea) (680402)
- IV pole (680401)
- Oxygen tank holder - D/E size cylinders (AL680403)
- 6 compartment drawer trays (3 ea)
(fit three per drawer), (40108)
- 7 compartment drawer tray (3 ea)
(fit three per drawer), (40109)

Drawer Configuration:

- (4) – 3" (7.6 cm)
- (1) – 6" (15.2 cm)
- (1) – 12" (30.5 cm) Storage compartment with flip up door

Internal Drawer Dimensions:

23" W x 17" D (58.4 cm W x 43.2 cm D)