

Classic Line, Nine Drawer, Pediatric Resuscitation Cart, Super-Stat Specialty Package #6411PEC



Construction:

- Sturdy, uni-body All Steel Construction

Features:

- Tall cabinet with separate breakaway locks on each drawer
- Pack of 100 break-away seals
- 30" of vertical drawer space
- Cart is designed to be compatible with popular pediatric color coding systems
- Includes multi-colored plates installed to identify drawers
- Stabilizing pontoon bumpers
- 5" ball-bearing casters, all full swivel - one tracking, one with brake
- Antimicrobial plastic top, removable
- Painted steel top rail
- Pre-drilled threaded holes for side mounted accessories
- Ball-bearing full extension drawer slides
- Latex free
- Durable, hard-baked powder coat finish
- Cart is delivered assembled, IV pole is customer installed
- Twelve year warranty on manufactured parts

Accessories Included:

- Cardiac board and brackets (680406)
- Three utility hooks
- IV pole (680401)
- Painted steel defibrillator shelf 20" x 12" (680323)
- Suction unit shelf (680305)
- Oxygen tank holder - D/E size cylinders (AL680403)
- Chart holder (680417)
- Locking sharps container, 5qt sharps, and glove box holder (684811)
- Two large utility baskets (680447)
- Hospital Grade electrical Outlet Strip with Cord Wrap not included as of January 15, 2015 due to UL specification change. Outlet conforming to UL1363A can be purchased as an accessory. (680404A)

Drawer Configuration:

- (8) – 3" (7.6 cm)
- (1) – 6" (15.2 cm)

Internal Drawer Dimensions:

- 23" W x 17" D (58.4 cm W x 43.2 cm D)

Overall Cart Dimensions:

- 45.75" H x 39" W x 24" D (116.2 cm H x 99 cm W x 60 cm D)
- With IV pole: 86.25" (219 cm)
- Writing surface: 41.25" (104.8cm)
- Cart weight: 187 pounds (84.8 kg)
- Approx. shipping weight: 241 lbs. (109.3 kg)
- Approx. shipping dims: 57" H x 33" W x 26" D (144.8cm H x 83.8cm W x 66cm D)