

## Classic Line, Six Drawer, Emergency Cart, Super Stat Accessory Package #6411



Color Shown is Red  
Additional Colors Available

### Drawer Configuration:

- (1) - 6-inch (15.2 cm)
- (3) - 3-inch (7.6 cm)
- (1) - 6-inch (15.2 cm)
- (1) - 9-inch (22.9 cm)

### Assembled Cart Dimensions

- 44.5" H x 47.5" W x 22" D  
(113.0 cm H x 120.7 W x 55.9 D)
- With IV pole: 85" (215.9 cm)
- Writing surface: 41.25" (104.8cm)
- Cart net weight: 191 lbs (86.8 kg)
- Approx. shipping weight: 245 lbs (111.4 kg)
- Approx. shipping dims: 57" H x 33" W x 26" D  
(144.8cm H x 83.8cm W x 66cm D)

### Construction:

- Sturdy, uni-body steel construction

### Features:

- Painted steel break-away lock bar
- Pack of 100 break-away seals
- 30" of vertical drawer space
- Cabinet and drawer design allows for customized drawer configurations
- Painted steel top rail / push handle
- Antimicrobial plastic top, removable
- Stabilizing pontoon bumpers
- Five-inch (12.7 cm) ball-bearing premium casters, all full swivel-one tracking, one braking
- Ball-bearing full extension drawer slides
- Cart is delivered mostly assembled
- Pre-drilled threaded holes for side mounted accessories
- Latex free
- Durable, hard-baked powder coat finish
- Twelve year warranty on manufactured parts

### Super Stat Accessory Package Includes:

- Cardiac board and brackets (680406)
- Three utility hooks (UHOOKS2)
- IV pole (680401)
- Painted steel defibrillator shelf 20" x 12" (680323)
- Suction unit shelf (SUCTIONSHLF)
- Oxygen tank holder - D/E size cylinders (O2HLDR)
- Chart holder (CHARTHLDRRC)
- Locking sharps container, 5qt sharps, and glove box holder (684811)
- Two large utility baskets (BASKETRC)
- Two 6 compartment drawer trays (fit three per drawer), (TRAY6COMP)
- One 7 compartment drawer trays (fit three per drawer), (TRAY7COMP)
- Two 3" drawer trays with adjustable dividers (fit two per drawer), (EXTRAY3)
- Four 6" drawer trays with adjustable dividers (fit two per drawer), (EXTRAY6)
- Hospital Grade electrical Outlet Strip with Cord Wrap not included as of January 15, 2015 due to UL specification change. Outlet conforming to UL1363A can be purchased as an accessory. (OUTLET6)