



650 Ford Street
 Colorado Springs, CO 80915
 Toll:(800)433-4064 * Ph: (719)637-0300
 Fax:(719)597-8273 * www.harloff.com

**Classic Line, Three Drawer Infection Control Cart with
 Breakaway Lock, Specialty Package
 #6501**



Specifications

Construction:

- Sturdy, uni-body steel construction

Features:

- Stainless steel break-away lock bar
- Pack of 100 break-away seals
- 30" of vertical drawer space
- Cabinet and drawer design allows for customized drawer configurations
- Full wrap-around bumper on separate stabilizing frame
- 5" ball-bearing casters, all full swivel - one tracking, one with brake
- Ball-bearing full extension drawer slides
- Removable plastic top
- Stainless steel top rail / push handle (68160)
- Coat and hat rack (680412)
- Utility hooks (4 ea) (680402)
- Cart is delivered mostly assembled
- Pre-drilled threaded holes for side mounted accessories
- Latex free
- Durable, hard-baked powder coat finish
- Twelve year warranty

Isolation Cart Upgrade:

- Silver ion antimicrobial paint available

Internal Drawer Dimensions:

- 23" W x 17" D (58.42 cm W x 43.18 cm D)

Assembled Cart Dimensions

- 44.5" H x 32" W x 22" D (113 cm H x 81.3 cm W x 55.9 cm D)
- With coat and hat rack: 62.5" H (158.8 cm H)
- Cart net weight: 120 pounds (63.6 kg)
- Approx. shipping weight: 150 pounds (72.7 kg)

Drawer Configuration:

- (1) - 6" (15.2 cm) multi-purpose storage drawer
- (2) - 12" (30.5 cm) multi-purpose storage drawer

Available in more than twenty color options:

View them online: <https://www.harloff.com/finish-options/>