

Emergency Anesthesia Cart/ Model: 6600

Classic Line, Six Drawer Emergency Anesthesia Cart with Breakaway Lock, Standard Package



Eighteen Available Colors:

Antimicrobial Finish Available: Lemon Drop Yellow Only
Beige, Burgundy, Light Blue, Navy Blue, Light Gray, Mauve, Red, Teal, Yellow, Forest Green, Hammertone Blue, Hammertone Gray, Hammertone Green, Hammer Red, Hammer Black, White, and Sand

Warranty:

- 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments.
- 1-Year Limited Warranty on any bought-out parts

Construction:

- Sturdy, uni-body steel construction

Features:

- Breakaway lock
- 30" of vertical drawer space
- Cabinet and drawer design allows for customized drawer configurations
- Ball-bearing full extension drawer slides
- 5" ball-bearing casters, all full swivel - one directional, one with brake
- Antimicrobial plastic top, removable
- Pre-drilled threaded holes for side mounted accessories
- Two aluminum side rails / push handles
- Stabilizing pontoon bumpers
- Latex free
- Durable, hard-baked powder coat finish

Drawer Configuration:

- (6) – 3" (7.6 cm)
- (2) – 6" (15.2 cm)

Internal Drawer Dimensions:

- 23" W x 17" D (58.4 cm W x 43.2 cm D)

Overall Cart Dimensions:

- 42.5" H x 32" W x 22" D (107.95 cm H x 81.3 W x 55.9 D)
- Writing surface: 41.25" H x (104.8 cm)
- Cart net weight: 129 lbs (58.5 kg)
- Approximate shipping weight: 170 lbs (77.1 kg)
- Approximate shipping dims: 57" H x 33" W x 26" D (144.8 cm H x 83.8 cm W x 66 cm D)