

OptimAI Line Cart Model: AL16CS-5

16 Bin Cassette Medication with Key Lock, Cassette Exchange Package



Seventeen Available Colors:

Antimicrobial Finish Available: Lemon Drop Yellow Only
Beige, Burgundy, Light Blue, Navy Blue, Light Gray, Mauve, Red, Teal, Yellow, Forest Green, Hammertone Blue, Hammertone Gray, Hammertone Green, Hammer Red, Hammer Black, White, and Sand

Warranty:

- 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments.
- 1-Year Limited Warranty on any bought-out parts

Construction:

- Painted aluminum cabinet and drawers

Features:

- 16 each 5" wide (12.7 cm) patient bin drawers in 2 removable cassettes; bin dividers and labels are included
- Cassettes are located in pull-out drawers for ergonomic removal
- BEST® lock key locking system with two keys
- Pull-out shelf on right side of cart
- Pre-threaded accessory mounting holes
- Self-closing, ball-bearing full extension drawer slides
- Removable plastic top with integrated push handles
- Corner bumpers, protecting the cart
- 5" ball-bearing casters, all full swivel - one directional, two with brake
- Latex Free
- Hard-baked powder coat finish

Suggested Cart Options:

- Cassette exchange package: 16CS5-PKG (includes 16 each 5" wide (12.7 cm) patient bin drawers and 2 double-tier aluminum cassettes)
- Narcotics drawer: CS-4FLN (separately locking narcotics drawer for 4" multi-purpose storage drawer)

Internal Drawer Dimensions:

- 20" W x 17" D (50.8 cm W x 43.2 cm D)

Bin Dimensions:

- 5" W x 9.5" L x 3" H (12.7 cm W x 24.1 cm L x 7.6 cm H)

Drawer Configuration:

- (1) - 4-inch (10.2 cm) multi-purpose storage drawer
- (2) - 8-inch (20.3 cm) double tier cassettes, with eight 5-inch wide patient bin drawers each
- (1) - 4-inch (20.3 cm) multi-purpose storage drawer
- (1) - 8-inch (20.3 cm) multi-purpose storage drawer

Overall Cart Dimensions:

- 45.7" H x 30.5" W x 23.8" D (116.2 cm H x 77.5 cm W x 60.6 cm D)
- Writing surface: 43.75" H (111.1 cm)
- Cart net weight: 132 lbs (59.9 kg)
- Approximate shipping weight: 180 lbs (81.8 kg)
- Approximate shipping dims: 48" H x 33" W x 26" D (121.9 cm H x 83.8 cm W x 66 cm D)