

## OptimAI Line Cart Model: AL24CS-5

24 Bin Cassette Medication with Key Lock, Cassette Exchange Package



### Seventeen Available Colors:

**Antimicrobial Finish Available: Lemon Drop Yellow Only**  
Beige, Burgundy, Light Blue, Navy Blue, Light Gray, Mauve, Red, Teal, Yellow, Forest Green, Hammertone Blue, Hammertone Gray, Hammertone Green, Hammer Red, Hammer Black, White, and Sand

### Warranty:

- 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments.
- 1-Year Limited Warranty on any bought-out parts

### Construction:

- Painted aluminum cabinet and drawers

### Features:

- 24 each 5" wide (12.7cm) patient bin drawers in 3 removable cassettes; bin dividers and labels are included
- Cassettes are located in pull-out drawers for ergonomic removal
- BEST® lock key locking system with two keys
- Pull-out shelf on right side of cart
- Pre-threaded accessory mounting holes
- Self-closing, ball-bearing full extension drawer slides
- Removable plastic top with integrated push handles
- Corner bumpers, protecting the cart
- 5" ball-bearing casters, all full swivel - one directional, two with brake
- Latex free
- Hard-baked powder coat finish

### Suggested Cart Options:

- Cassette exchange package: 24CS5-PKG (includes 24 each 5" wide (12.7 cm) patient bin drawers and 3 double-tier aluminum cassettes)
- Narcotics drawer: CS-4FLN (separately locking narcotics drawer for 4" multi-purpose storage drawer)

### Internal Drawer Dimensions:

- 20" W x 17" D (50.8 cm W x 43.2 cm D)

### Bin Dimensions:

- 5" W x 9.5" L x 3" H (12.7 cm W x 24.1 cm L x 7.6 cm H)

### Drawer Configuration:

- (1) - 4-inch (10.2 cm) multi-purpose storage drawer
- (3) - 8-inch (20.3 cm) double tier cassettes, with twelve 3.5-inch wide patient bin drawers each
- (1) - 4-inch (20.3 cm) multi-purpose storage drawer

### Overall Cart Dimensions:

- 45.7" H x 30.5" W x 23.8" D (116.2 cm H x 77.5 cm W x 60.6 cm D)
- Writing surface: 43.75" H (111.1 cm)
- Cart net weight: 132 lbs (59.9 kg)
- Approximate shipping weight: 180 lbs (81.8 kg)
- Approximate shipping dims: 48" H x 33" W x 26" D (121.9 cm H x 83.8 cm W x 66 cm D)