

Aluminum Mini Line Short Cart/ Model: AL3233E

Three Drawer Aluminum with Basic Electronic Pushbutton Lock, Standard Package



Seventeen Available Colors:

Antimicrobial Finish Available: Lemon Drop Yellow Only
Beige, Burgundy, Light Blue, Navy Blue, Light Gray, Mauve, Red, Teal, Yellow, Forest Green, Hammertone Blue, Hammertone Gray, Hammertone Green, Hammer Red, Hammer Black, White, and Sand

Construction:

- Painted Aluminum cabinet and drawers

Features:

- Space-saver design fits into smaller areas
- 5-inch ball bearing casters, full swivel, 2 with brake
- Basic electronic pushbutton lock
- Key lock override
- Factory set user and supervisor code
- Cabinet may be programmed in the field
- Stabilizing pontoon bumper
- Pre-threaded accessory mounting holes
- Ball-bearing full extension drawer slides
- Antimicrobial plastic top, removable
- Latex free
- Durable hard-baked powder coat finish

Drawer Configuration:

- (1) - 3-inch (7.6 cm)
- (1) - 6-inch (15.2 cm)
- (1) - 9-inch (22.9 cm)

Internal Drawer Dimensions:

- 12.6" W x 16.75" D (38.1 cm W x 42.5 cm D)

Overall Cart Dimensions:

- 21.50" W x 22.50" D x 30" H (54.61 cm W x 57.15 cm D x 76.20 cm H)
- Writing surface: 29" H (73.7 cm)
- Cart net weight: 76 lbs (34.4 kg)
- Approximate shipping weight: 96 lbs (43.5 kg)
- Approximate shipping dims: 24"W x 29.5"D x 42.5"H (60.9 cm W x 74.9 cm D x 107.9 cm H)

Warranty:

- 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments.
- 1-Year Limited Warranty on any bought-out parts