

Aluminum Mini Line Short Cart/ Model: AL3233ELP

Three Drawer Aluminum with Electronic Lock and Proximity Card and Keypad Access, Standard Package



Seventeen Available Colors:

Antimicrobial Finish Available: Lemon Drop Yellow Only
Beige, Burgundy, Light Blue, Navy Blue, Light Gray, Mauve, Red, Teal, Yellow, Forest Green, Hammertone Blue, Hammertone Gray, Hammertone Green, Hammer Red, Hammer Black, White, and Sand

Warranty:

- 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments.
- 1-Year Limited Warranty on any bought-out parts

Construction:

- Painted Aluminum cabinet and drawers

Features:

- Space-saver design fits into smaller areas
- 5-inch ball bearing casters, full swivel, 2 with brake
- Electronic-locking cart with proximity card, keypad access, and key lock override with 2 keys
- Drawers automatically re-lock after access (pre-set or programmed time)
- Access rights may be assigned for up to 250 users per cart
- Allows for supervisor assignments; dual credential access available
- May add, edit, or delete users and user groups
- Battery-operated with six D-cell batteries; life of one year (minimum) at 25 to 50 uses/day
- Lock has an LED low battery indicator
- Stabilizing pontoon bumper
- Pre-threaded accessory mounting holes
- Ball-bearing full extension drawer slides
- Antimicrobial plastic top, removable
- Latex free
- Durable hard-baked powder coat finish

Drawer Configuration:

- (1) - 3-inch (7.6 cm)
- (1) - 6-inch (15.2 cm)
- (1) - 9-inch (22.9 cm)

Internal Drawer Dimensions:

- 12.6" W x 16.75" D (38.1 cm W x 42.5 cm D)

Overall Cart Dimensions:

- 21.50" W x 22.50" D x 30" H (54.61 cm W x 57.15 cm D x 76.20 cm H)
- Writing surface: 29" H (73.7 cm)
- Cart net weight: 76 lbs (34.4 kg)
- Approximate shipping weight: 96 lbs (43.5 kg)
- Approximate shipping dims: 24"W x 29.5"D x 42.5"H (60.9 cm W x 74.9 cm D x 107.9 cm H)