

Aluminum Mini Line Short Cart/ Model: AL3235EKC

Five Drawer Aluminum with CompX® Electronic Keypad Lock, Standard Package



Seventeen Available Colors:

Antimicrobial Finish Available: Lemon Drop Yellow Only
Beige, Burgundy, Light Blue, Navy Blue, Light Gray, Mauve, Red, Teal, Yellow, Forest Green, Hammertone Blue, Hammertone Gray, Hammertone Green, Hammer Red, Hammer Black, White, and Sand

Construction:

- Painted Aluminum cabinet and drawers

Features:

- Space-saver design fits into smaller areas
- 5-inch ball bearing casters, full swivel, 2 with brake
- Electronic locking cart with keypad access and Key Lock override
- Drawers automatically re-lock after access (pre-set or programmed time)
- Access rights may be assigned for up to 250 users per cart
- May add, edit, or delete user and user groups
- Battery operated with six D-cell Batteries
- Lock has an LED low battery indicator
- Stabilizing pontoon bumper
- Pre-threaded accessory mounting holes
- Ball-bearing full extension drawer slides
- Antimicrobial plastic top, removable
- Latex free
- Durable hard-baked powder coat finish

Drawer Configuration:

- (4) - 3-inch (7.6 cm)
- (1) - 6-inch (15.2 cm)

Internal Drawer Dimensions:

- 12.6" W x 16.75" D (38.1 cm W x 42.5 cm D)

Overall Cart Dimensions:

- 21.50" W x 22.50" D x 30" H (54.61 cm W x 57.15 cm D x 76.20 cm H)
- Writing surface: 29" H (73.7 cm)
- Cart net weight: 76 lbs (34.4 kg)
- Approximate shipping weight: 96 lbs (43.5 kg)
- Approximate shipping dims: 24"W x 29.5"D x 42.5"H (60.9 cm W x 74.9 cm D x 107.9 cm H)

Warranty:

- 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments.
- 1-Year Limited Warranty on any bought-out parts