

OptimAL Line 600/ Model: AL602PC

Aluminum Punch Card Medication Cart, Standard Package, Double Wide Narcotics Drawer

OptimAL Line medication carts are lightweight due to the all-aluminum construction. These carts are available to accommodate all types of medication packaging and offer a complete line of accessories to meet your storage and security needs. Harloff carts are shipped fully assembled and ready to use.



Seventeen Available Colors:

Antimicrobial Finish Available: Lemon Drop Yellow Only
Beige, Burgundy, Light Blue, Navy Blue, Light Gray, Mauve, Red, Teal, Yellow, Forest Green, Hammertone Blue, Hammertone Gray, Hammertone Green, Hammer Red, Hammer Black, White, and Sand

Construction:

- Painted aluminum cabinet and drawers

Features:

- Dual column medication cart with 600 punch card capacity, double width narcotics drawer
- BEST® brand lock on cabinet, 4 keys
- Illinois lock on narcotics drawer, 4 keys
- Painted pull-out shelf
- Ball-bearing, full extension drawer slides
- Plastic replaceable top with integrated push handles
- Corner bumpers, protecting the cart
- 5" ball-bearing premium rubber casters, full swivel, 2 with toe locking brake and 1 directional
- Hard-baked powder coat finish

Left Column Drawer Configuration:

- (1) - 4-inch multipurpose storage drawer
- (3) - 9.75-inch punch card drawers

Right Column Drawer Configuration:

- (1) - 4-inch multipurpose storage drawer
- (1) - 9.75-inch locking narcotics drawer
- (2) - 9.75-inch multipurpose storage drawer

Assembled Cart Dimensions:

- 45.75" H x 43.75" W x 24" D (116.2 cm H x 109.2 cm W x 59.7 cm D)
- Writing surface: 44" (111.8 cm)
- Cart weight: 145 pounds (65.9 kg)
- Approx. shipping weight: 190 pounds (86.36 kg)
- Approx. shipping dims: 57" H x 47" W x 24" D (144.8 cm H x 119.4 cm W x 60.9 cm D)

Warranty:

- 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments.
- 1-Year Limited Warranty on any bought-out parts