



650 Ford Street
Colorado Springs, CO 80915
Toll:(800)433-4064 * Ph: (719)637-0300
Fax:(719)597-8273 * www.harloff.com

Aluminum Universal Line Tall 6 Drawer Crash Cart with Breakaway Lock and Emergency Accessory Package #AL809B6-EMG

Harloff's 20 years of innovation in hospital cart design and construction has resulted in the Universal Line. Its modern appearance incorporates industry-leading standard features with flexible configuration options.



Specifications

Construction:

- Aluminum construction

Features:

- Internal pullout shelf on left side
- Break-away lock bar and 100 numbered break-away seals
- 29.25" of vertical drawer space
- Cabinet design allows for customized drawer configurations
- Painted steel top rail / push handle
- Removable plastic top
- Protective pontoon bumpers
- Extended caster footprint for stability when drawers are extended
- 5" (12.7 cm) ball-bearing premium casters, all full swivel
 - one tracking, one braking
- Ball-bearing full extension drawer slides
- Cart is delivered wheels-on and mostly assembled
- Pre-drilled hole locations for standard accessories – many are pre-threaded for ease of installation
- Latex free
- Durable, hard-baked powder coat finish available in 16 colors
- Twelve year manufacturer's warranty

Package Accessories Include:

- Cardiac board and brackets
- Utility hooks (2 ea)
- IV pole
- Aluminum oxygen tank holder - D/E size cylinders
- Full Drawer Trays with adjustable divider set, 68530-P1 (2 ea)

Drawer Configuration:

- (4) – 3.25" (8.3 cm)
- (1) – 6.5" (16.5 cm)
- (1) – 9.75" (24.8 cm)

Internal Drawer Dimensions:

- 23" W x 17" D (58.4 cm W x 43.2 cm D)

Sixteen available color choices:

*Antimicrobial Finish Available

Burgundy, Hammer Red, Red, Yellow, Beige, Sand, White, HammerTone Green, Teal, HammerTone Blue, Navy, Light Blue, HammerTone Gray, Light Gray, Mauve, Hammer Black

Assembled Cart Dimensions:

- 67.75" H x 35" W x 24.7" D (172.1 cm H x 88.9 cm W x 62.9 cm D)
- IV pole extends to 85.25" H (216 cm)
- Cart weight: 124 pounds (56.2 kg)
- Approximate shipping weight: 142 pounds (64.4 kg)