

Procedure Cart/ Model: AL809K6

Aluminum Universal Line Tall 6 Drawer Procedure Cart with Key Lock

Harloff's innovation in hospital cart design and construction has resulted in the Universal Line. Its modern appearance incorporates industry-leading standard features with flexible configuration options.



Seventeen Available Colors:

Antimicrobial Finish Available: Lemon Drop Yellow Only
Beige, Burgundy, Light Blue, Navy Blue, Light Gray, Mauve, Red, Teal, Yellow, Forest Green, Hammertone Blue, Hammertone Gray, Hammertone Green, Hammer Red, Hammer Black, White, and Sand

Construction:

- Painted aluminum cabinet and drawers

Features:

- Internal pullout shelf on right side
- Key lock with two keys
- 29.25" of vertical drawer space
- Cabinet design allows for customized drawer configurations
- Painted aluminum top rail / push handle
- Antimicrobial plastic top, removable
- Stabilizing pontoon bumpers
- 5-inch ball bearing casters, full swivel, 2 with brake
- Ball-bearing full extension drawer slides
- Pre-drilled threaded holes for standard accessories
- Latex free
- Durable hard-baked powder coat finish

Drawer Configuration:

- (4) - 3.25-inch (8.26 cm)
- (1) - 6.5-inch (16.5 cm)
- (1) - 9.75-inch (25.4 cm)

Internal Drawer Dimensions:

- 23" W x 17" D (58.4 cm W x 43.2 cm D)

Overall Cabinet Dimensions:

- 44.5" H x 31.5" W x 24.75" D (113 cm H x 80.8 cm W x 62.9 cm D)
- Writing surface: 43.5" H x (110.4 cm)
- Cart net weight: 104 lbs (47.3 kg)
- Approximate shipping weight: 130 lbs (59.1 kg)
- Approximate shipping dims: 56" H x 36" W x 26" D (142.2 cm H x 91.4 cm W x 66 cm D)

Warranty:

- 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments.
- 1-Year Limited Warranty on any bought-out parts