

## Universal Line Aluminum Super 6 Drawer Cart with Breakaway Lock and Super-Stat Emergency Package #AL810B6-EMG3

Harloff's 20 years of innovation in hospital cart design and construction has resulted in the Universal Line. Its modern appearance incorporates industry-leading standard features with flexible configuration options.



Color Shown is RED.  
Additional Colors Available

### Assembled Cart Dimensions:

- 52" H x 45.5" W x 25" D (132 cm H x 115.6 cm W x 63.5 cm D)
- IV Pole Extends to: 88.5" (225 cm) H
- Cart weight: 145 pounds (66 kg)
- Approximate shipping weight: 188 pounds (85.5 kg)

### Construction:

- All aluminum, light-weight construction

### Features:

- Internal Pullout shelf on right side side
- Break-away lock bar and 100 numbered break-away seals
- 32.5" of vertical drawer space
- Cabinet design allows for customized drawer configurations
- Painted steel top rail / push handle
- Removable plastic top
- Extended caster footprint for stability with protective bumper
- 5" (12.7 cm) ball-bearing premium casters, all full swivel-one tracking, one braking
- Ball-bearing full extension drawer slides
- Cart is delivered wheels-on and mostly assembled
- Pre-drilled hole locations for standard accessories – many are pre-threaded for ease of installation
- Latex free
- Durable, hard-baked powder coat finish available in 17 colors
- Twelve year manufacturer's warranty

### Super Stat Emergency Package Includes:

- Cardiac board and brackets (UL-CARDBRD)
- Three utility hooks
- IV pole (UL-IV2)
- Painted steel defibrillator shelf 20" x 12" (UL-DEFIB)
- Suction unit shelf
- Oxygen tank holder - D/E size cylinders (UL-ALO2)
- Chart holder (680416)
- Locking sharps container, 5qt sharps, and glove box holder (684811)
- Two large utility baskets (681501)
- Two 6 compartment drawer trays (fit three per drawer), (40108)
- One 7 compartment drawer trays (fit three per drawer), (40109)
- Two 3" drawer trays with adjustable dividers (fit two per drawer), (681501)
- Four 6" drawer trays with adjustable dividers (fit two per drawer), (681502)

• Hospital Grade Electrical Outlet Strip with Cord Wrap not included as of January 15, 2015 due to UL Specification change. Outlets conforming to UL1363A can be purchased as an accessory item. (Part #UL-OUTLETA)

### Internal Drawer Dimensions:

- 23" W x 17" D (58.4 cm W x 43.2 cm D)

### Drawer Configuration:

- (3) – 3.25" (8.3 cm)
- (2) – 6.5" (16.5 cm)
- (1) – 9.75" (24.8 cm)