



MED CARTS

AL8W10EKC4DSP

DISPILL® MEDICATION CART

132 Card Capacity & Electronic Lock w/Keypad

Features:

- Single column medication delivery cart for Dispill® cards and liquid storage with a maximum capacity of 132 Dispill® cards
- Includes one single width locking narcotics box made for Dispill® packaging that secures with a standard key lock
- Electronic locking cart with keypad access and key lock override
- Drawers automatically re-lock after access (pre-set or programmed time)
- Access rights may be assigned for up to 250 users per cart
- May add, edit, or delete user & user groups
- Uses six D-cell batteries; 1 year minimum life at 25-50 uses per day
- Lock has an LED low battery indicator
- All-aluminum construction with scratch resistant heat sealed powder coat finish
- Pull-out shelf
- Ball-bearing, full extension drawer slides
- Removeable plastic top
- 5" ball-bearing premium rubber casters, full swivel, two with toe locking brake and one directional
- Painted aluminum top rail / push handle
- Protective pontoon bumpers
- Extended caster footprint for stability when drawers are extended
- Latex free cart is delivered fully assembled
- Body has pre-threaded accessory mounting holes



Color shown is Hammertone Green. More than 20 additional colors are available.

Aluminum 132 Dispill® Card Capacity Medication Cart with Key Lock on Cabinet

Designed to meet the storage needs of an individual, a group home or a large nursing unit. Our steel, aluminum, wood look and value lines offer medication storage solutions at price points to meet your budget.

Internal Drawer Dimensions:

- 30" W x 17" D

Assembled Cart Dimensions:

- Writing Surface: 44"
- Cart weight: 130 pounds (59 kg)
- 48.5" H x 38.25" W x 25" D (123.2 cm x 97.2 cm x 63.5 cm)
- Approx shipping weight: 170 pounds (77.1 kg)
- Approx shipping dims: 52" H x 46" W x 26" D (132.1 cm x 116.8 cm x 66 cm)

Drawer Configuration:

- 2 ea – 3.25" multipurpose storage drawer
- 1 ea – 13" drawer with locking narcotics box & rail system to hold Dispill® packs
- 1 ea – 13" drawer with rail system to hold Dispill® packs

Warranty: 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts



Engineered & Manufactured in the U.S.



View similar products online at <https://www.harloff.com/product-category/specialty-packaging-carts/>

Toll Free: (800) 433-4064 • medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300