

ETC Line, Four Drawer, Treatment Cart #ETC-4



Over 20 Available Colors:

Click or copy & paste the link below to view all current available color options:

<https://www.harloff.com/finish-options/>

Specifications

Construction:

- Sturdy, uni-body steel construction

Features:

- Low height make these carts ideal for bedside or under a counter
- Key lock with two keys
- 3" stem casters, two with brake
- Roller drawer slides
- Removable plastic top
- Hard-baked powder coat finish
- Certified latex free
- Twelve year warranty

Drawer Configuration:

- (1) – 2.75" (6.98 cm)
- (3) – 5.75" (14.6 cm)

Internal Drawer Dimensions:

- 20" W x 12.75" D (50.8cm W x 32.4 cm D)

Assembled Cart Dimensions

- 28.63" H x 23.5" W x 15" D (72.73 cm H x 59.4 W x 38.4 D)
- Cart weight: 55 lbs (25 kg)
- Approx. Shipping Weight: 86 pounds (39.1 kg)
- Approx. Shipping Dims: 33" H x 29" W x 23" D (83.8 cm H x 73.7 cm W x 58.4 cm D)

** PLEASE NOTE THIS CART IS SHORTER THAN STANDARD COUNTER HEIGHTS**