

## IntraOcular Lens Storage Cart/ Model: IOL1260

### Six Drawers with Key Lock

Attractive, secure storage for IntraOcular lens. Provides mobility and is compatible with multiple brands and sizes of IOL boxes, including Alcon, Bausch and Lomb, Rayner, and Abbot Medical Optics.



### Construction:

- Painted steel cabinet and drawers

### Features:

- Lockable, secure storage of supplies: key lock with two keys
- Contemporary design with rounded corners
- Integrated side push handles
- Six pull-out drawer shelves with adjustable dividers – ideal for medical supplies used in surgery centers and operating rooms
- Double doors with acrylic inserts open 270 degrees for full access
- Acrylic windows allow for visibility of supplies, while locking doors provide security
- Ball-bearing full extension drawer slides
- 5" ball-bearing casters, all full swivel - one directional, one with brake
- Removable plastic top with integrated push handles
- Protective corner bumpers
- Latex free

### Internal Drawer Dimensions:

- 33.75" W x 24" D x 3.63" H  
(85.5 cm x 61 cm x 9.22 cm)
- 7" Vertical spacing between drawers

### Overall Cart Dimensions:

- 55.5" H x 43.4" W x 30.75" D  
(141 cm H x 110.2 W x 78.1 D)
- Cart weight: 385 lbs (175 kg)
- Approximate shipping weight: 430 lbs (209 kg)
- Approximate shipping dims: 62" H x 48" W x 33.25" D  
(157.5 cm H x 121.9 cm W x 84.5" cm W)

### Colors:

More than twenty color options are available.

### Choose Your Colors Online:

<https://www.harloff.com/finish-options/>

### Warranty:

- 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments.
- 1-Year Limited Warranty on any bought-out parts