

IntraOcular Lens Storage Cart/ Model: IOL924

Six Drawers with Key Lock

Attractive, secure storage for IntraOcular lens. Provides mobility and is compatible with multiple brands and sizes of IOL boxes, including Alcon, Bausch and Lomb, Rayner, and Abbot Medical Optics.



Seventeen Available Colors:

Antimicrobial Finish Available: Lemon Drop Yellow Only
Beige, Burgundy, Light Blue, Navy Blue, Light Gray, Mauve, Red, Teal, Yellow, Forest Green, Hammertone Blue, Hammertone Gray, Hammertone Green, Hammer Red, Hammer Black, White, and Sand

Construction:

- Painted steel cabinet and drawers

Features:

- Lockable, secure storage of supplies: key lock with two keys
- Contemporary design with rounded corners
- Six pull-out drawer shelves with adjustable dividers – ideal for medical supplies used in surgery centers and operating rooms
- Double doors with plexi-glass inserts open 180 degrees for full access
- Plexi-glass windows allow for visibility of supplies, while locking doors provide security
- Ball-bearing full extension drawer slides
- 5" ball-bearing casters, all full swivel - one directional, one with brake
- Removable plastic top with integrated push handles
- Protective corner bumpers
- Latex free

Internal Drawer Dimensions:

- 33.75" W x 17.8" D x 3.63" H
(85.5 cm x 45.2 cm x 9.22 cm)
- 7" Vertical spacing between drawers

Overall Cart Dimensions:

- 55.5" H x 43.4" W x 23.5" D
(141 cm H x 110.2 W x 60.7 D)
- Cart weight: 277 lbs (175 kg)
- Approximate shipping weight: 334 lbs (151.8 kg)
- Approximate shipping dims: 62" H x 48" W x 33.25" D
(157.5 cm H x 121.9 cm W x 84.5" cm W)

Warranty:

- 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments.
- 1-Year Limited Warranty on any bought-out parts