



Wood Style

M3DS1821KC4-AN

M-SERIES X-SHORT WOOD STYLE CART

with Four Drawers and a Standard Key Lock

Why Wood Style Drawer Fronts?

Medical carts which appear more like furniture tend to provide better psychological benefits for your patients, as they feel more at home and comfortable than they would in a sterile setting.

Features:

- **Wood style drawer fronts giving the cart a furniture feel are available in multiple finish options**
- Sturdy, uni-body steel construction
- Key lock with two keys
- Durable scratch resistant heat sealed powder coat finish on body
- Full extension drawers slides for easy access
- Pre-threaded mounting holes for accessories found at harloff.com/product-category/accessories/
- 21" of vertical drawer space
- Cart is delivered assembled
- Antimicrobial plastic top, removable
- Fully customizable drawer configurations
- 3" (7.6 cm) single wheel stem caster set, full swivel - two with brake

Two standard wood vinyl colors available plus additional special order vinyl color options with a minimum cart purchase and extended lead time.



Colors:

Body shown in White, but available in more than twenty color options which you can view online at www.harloff.com/finish-options/. Drawers shown in Wood Vinyl American Natural color.



Engineered & Manufactured in the U.S.

Overall Cart Dimensions:

- Internal Drawer Dimensions: 12.5" W x 17" D (31.8 cm W x 43.2 cm D)
- Work Surface Height: 30.25" (76.8 cm)
- Cart dims: 31.375" H x 18" W x 18" D (79.7 cm H x 45.7 cm W x 45.7 cm D)
- Cart net weight: 78 lbs (35.4 kg)
- Approx shipping weight: 100 lbs (45.4 kg)
- Approx shipping dims: 41" H x 29" W x 24" D (104 cm H x 73.7 cm W x 61 cm D)

Drawer Configuration:

- One 3" drawers (7.6 cm) and three 6" drawers (15.2 cm)

Warranty: 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts

Hospitals, Surgery Centers and Clinics can bring the warm appearance of wood to their treatment and other procedure carts with Harloff's line of M-Series carts with vinyl-coated aluminum drawers. Designed and built to combine the aesthetics of wood with the long-term functionality of a Harloff metal cart.

View similar products online at

<https://www.harloff.com/product-category/m-series-medical-carts/>

**Toll Free: (800) 433-4064 • medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300**

HARLOFF
INNOVATING STORAGE FOR CARE