



M-series

M3DS1824B03

M-SERIES SHORT MINI CART

Three Drawers with a Breakaway Lock

Features:

- Sturdy, uni-body steel construction
- Heat sealed powder coat steel break-away lock bar
- Pack of 100 breakaway seals
- Durable scratch resistant heat sealed powder coat finish
- Full extension drawers slides for easy access
- Body has pre-threaded accessory mounting holes
- Antimicrobial plastic top, removable
- 24" of vertical drawer space
- Cart is delivered assembled
- Fully customizable drawer configurations
- 3" (7.6 cm) single wheel stem caster set, full swivel - two with brake
- Latex free



Color shown in Red. More than 20 additional colors are available.

Overall Cart Dimensions:

- Internal Drawer Dimensions: 12.5" W x 17" D (31.8 cm W x 43.2 cm D)
- Work Surface Height: 33.375" (84.8 cm)
- Cart dims: 34.5" H x 18" W x 18" D (87.6 cm H x 45.7 cm W x 45.7 cm D)
- Cart net weight: 83 lbs (37.6 kg)
- Approx shipping weight: 107 lbs (48.5 kg)
- Approx shipping dims: 50" H x 24" W x 24" D (127 cm H x 61 cm W x 61 cm D)

Drawer Configuration:

One 6" drawer (15.2 cm) and two 9" drawers (22.8 cm)

Warranty: 12-Year Manufacturer's Warranty on shipments to US and Canada. 1-Year Limited Warranty on any bought-out parts.

Short Height, Mini Width, Three Drawer Crash Cart with a Breakaway Lock

Designed and built to meet the needs of today's healthcare market. Comfort drawer pulls, expanded lock options and a standard antimicrobial top make things easier for your staff and help support your infection control initiatives. M-Series carts offer you choices.



Engineered & Manufactured in the U.S.

HARLOFF

INNOVATING STORAGE FOR CARE

View similar products online at

<https://www.harloff.com/product-category/emergency-crash-carts/>

Toll Free: (800) 433-4064 • medicalsales@harloff.com

www.harloff.com • Ph: (719) 637-0300