



# M-series

## M3DS1824E04

### M-SERIES SHORT MINI CART

Four Drawers with an Electronic Pushbutton Lock

#### Features:

- Sturdy, uni-body steel construction
- Basic electronic pushbutton lock
- Key lock override
- Factory set user and supervisor code
- Cabinet may be programmed in the field
- Durable scratch resistant heat sealed powder coat finish
- Full extension drawers slides for easy access
- Body has pre-threaded accessory mounting holes
- Antimicrobial plastic top, removable
- 24" of vertical drawer space
- Cart is delivered assembled
- Fully customizable drawer configurations
- 3" (7.6 cm) single wheel stem caster set, full swivel - two with brake
- Latex free

#### Overall Cart Dimensions:

- Internal Drawer Dimensions: 12.5" W x 17" D (31.8 cm W x 43.2 cm D)
- Work Surface Height: 33.375" (84.8 cm)
- Cart dims: 34.5" H x 18" W x 18" D (87.6 cm H x 45.7 cm W x 45.7 cm D)
- Cart net weight: 83 lbs (37.6 kg)
- Approx shipping weight: 107 lbs (48.5 kg)
- Approx shipping dims: 50" H x 24" W x 24" D (127 cm H x 61 cm W x 61 cm D)

#### Drawer Configuration:

One 3" drawer (7.6 cm), two 6" drawers (15.2 cm) and one 9" drawer (22.8 cm)

**Warranty:** 12-Year Manufacturer's Warranty on shipments to US and Canada. 1-Year Limited Warranty on any bought-out parts.



Color shown in Yellow. More than 20 additional colors are available.



Engineered & Manufactured in the U.S.

### Short Height, Mini Width, Four Drawer Infection Control Cart with Electronic Pushbutton Lock

Designed and built to meet the needs of today's healthcare market. Comfort drawer pulls, expanded lock options and a standard antimicrobial top make things easier for your staff and help support your infection control initiatives. M-Series carts offer you choices.



View similar products online at <https://www.harloff.com/product-category/m-series-medical-carts/>

Toll Free: (800) 433-4064 • [medicalsales@harloff.com](mailto:medicalsales@harloff.com)  
[www.harloff.com](http://www.harloff.com) • Ph: (719) 637-0300