



M-series

M3DS1824K05+MD18-ANS

M-SERIES SHORT ANESTHESIA CART

Mini 5-Drawer Cart with MD18-ANS Accessory Package

Short Height, Mini Width, Five Drawer Medical Cart with a Standard Key Lock

Designed and built to meet the needs of today's healthcare market. Ergonomic side handles, comfort drawer pulls, expanded lock options and a standard antimicrobial top make things easier for your staff and help support your infection control initiatives. M-Series carts offer you choices.

MD18-ANS Accessories Include:



3-Tier Back Rail
(MD18-RBR3)



(3) 12" Side Rails
(SIDERAIL12)



5-Compartment
Tilt Bin w/Rail Clips
(TILTBIN5RC)



Engineered &
Manufactured
in the U.S.



Features:

- Sturdy, uni-body steel construction
- Key lock with two keys
- Durable scratch resistant heat sealed powder coat finish
- Full extension drawers slides for easy access
- Body has pre-threaded accessory mounting holes
- Antimicrobial plastic top, removable
- 24" of vertical drawer space
- Cart is delivered assembled
- Fully customizable drawer configurations
- 3" (7.6 cm) single wheel stem caster set, full swivel - two with brake
- Latex free

Warranty:

- 12-Year Manufacturer's Warranty on shipments to US and Canada. 1-Year Limited Warranty on any bought-out parts.

HARLOFF
INNOVATING STORAGE FOR CARE

View similar products online at

<https://www.harloff.com/product-category/anesthesia-carts/>

Toll Free: (800) 433-4064 • Ph: (719) 637-0300
medicalsales@harloff.com • www.harloff.com

Drawer Configuration:

- 12.5" W x 17" D (31.8 cm W x 43.2 cm D)
- Three 3" drawers (7.6 cm), one 6" drawer (15.2 cm) and one 9" drawer (22.8 cm)

Assembled Cart Dimensions (without accessories):

- 34.5" H x 18" W x 18" D (87.6 cm H x 45.7 cm W x 45.7 cm D)
- Cart weight: 83 lbs (37.6 kg)
- Approximate shipping weight: 107 lbs (48.5 kg)
- Approximate Shipping Dims: 50" H x 24" W x 24" D (127 cm x 61 cm x 61 cm)