



M-series

M3DS1830E05

M-SERIES TALL MINI CART

Five Drawers with an Electronic Pushbutton Lock

Features:

- Sturdy, uni-body steel construction
- Basic electronic pushbutton lock
- Key lock override
- Factory set user and supervisor code
- Cabinet may be programmed in the field
- Durable scratch resistant heat sealed powder coat finish
- Full extension drawers slides for easy access
- Body has pre-threaded accessory mounting holes
- Antimicrobial plastic top, removable
- 30" of vertical drawer space
- Cart is delivered assembled
- Fully customizable drawer configurations
- 3" (7.6 cm) single wheel stem caster set, full swivel - two with brake
- Latex free



Engineered & Manufactured in the U.S.



Colors:

Shown in Light Blue, but available in more than twenty color options which you can view online at www.harloff.com/finish-options/.

Overall Cart Dimensions:

- Internal Drawer Dimensions: 12.5" W x 17" D (31.8 cm W x 43.2 cm D)
- Work Surface Height: 39.625" (100.6 cm)
- Cart dims: 40.75" H x 18" W x 18" D (103.5 cm H x 45.7 cm W x 45.7 cm D)
- Cart net weight: 95 lbs (43.1 kg)
- Approx shipping weight: 120 lbs (54.4 kg)
- Approx shipping dims: 50" H x 29" W x 24" D (127 cm H x 73.7 cm W x 61 cm D)

Drawer Configuration:

- Two 3" drawers (7.6 cm), one 6" drawer (15.2 cm) and two 9" drawers (22.8 cm)

Warranty: 12-Year Manufacturer's Warranty on shipments to US and Canada, 1-Year Limited Warranty on any bought-out parts

Tall Height, Mini Width, Five Drawer Procedure Cart with an Electronic Pushbutton Lock

Designed and built to meet the needs of today's healthcare market. Comfort drawer pulls, expanded lock options and a standard antimicrobial top make things easier for your staff and help support your infection control initiatives. M-Series carts offer you choices.



View similar products online at

<https://www.harloff.com/product-category/m-series-medical-carts/>

Toll Free: (800) 433-4064 • medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300