



**M-series**

**M3DS1830K06+MD18-ANS**

## M-SERIES TALL ANESTHESIA CART

### Mini Six Drawer Cart w/MD18-ANS Accessory Package

### Tall Height, Mini Width, Six Drawer Anesthesia Cart with a Standard Key Lock

Designed and built to meet the needs of today's healthcare market. Comfort drawer pulls, expanded lock options and a standard plastic top make things easier for your staff and help support your infection control initiatives. M-Series carts offer you choices.

### MD18-ANS Accessories Include:



3-Tier Back Rail  
(MD18-RBR3)



(3) 12" Side Rails  
(SIDERAIL12)



5-Compartment  
Tilt Bin w/Rail Clips  
(TILTBIN5RC)



Engineered &  
Manufactured  
in the U.S.



### Features:

- Sturdy, uni-body steel construction
- Key lock with two keys
- Durable scratch resistant heat sealed powder coat finish
- Full extension drawers slides for easy access
- Body has pre-threaded accessory mounting holes
- Removable plastic top
- 30" of vertical drawer space
- Drawer Configuration: Two 3" drawers (7.6 cm), one 6" drawer (15.2 cm) and two 9" drawers (22.8 cm)
- Cart is delivered assembled
- Fully customizable drawer configurations
- 3" (7.6 cm) single wheel stem caster set, full swivel - two with brake
- 12-Year Manufacturer's Warranty on shipments to US and Canada. 1-Year Limited Warranty on any bought-out parts.

**HARLOFF**  
INNOVATING STORAGE FOR CARE

View similar products online at

<https://www.harloff.com/product-category/anesthesia-carts/>

**Toll Free: (800) 433-4064 • Ph: (719) 637-0300**  
**medicalsales@harloff.com • www.harloff.com**

Color shown with a Cream body  
and Navy drawers. More than 20  
additional colors are available.

All dimensions are available in the "MA-Dims-18wide-3DScaster.pdf" document. Shipping weights can be found at the online product listing for this model, or feel free to request them from your Harloff Team.