



M-series M3DS2424B05+MD24-EMG

M-SERIES SHORT EMERGENCY CART with MD24-EMG Emergency Package

Short, Medium Width, Five Drawer Emergency Cart with Breakaway Lock

Designed and built to meet the needs of today's healthcare market. Ergonomic side handles, comfort drawer pulls, expanded lock options and a standard plastic top make things easier for your staff and help support your infection control initiatives. M-Series carts offer you choices.

MD24-EMG Accessories Include:



Oxygen tank holder (O2HLDR)



IV pole (MD-IV2)



Cardiac board and brackets (MD-CARDBRD)

Color shown is Red. More than 20 additional colors are available.



Engineered & Manufactured in the U.S.

Features:

- Sturdy, uni-body steel construction
- Choice of lock options
- Heat sealed powder coat steel breakaway lock bar
- Pack of 100 breakaway seals
- 24" of vertical drawer space
- Cart is delivered assembled
- Removable plastic top
- Latex free
- 3" (7.6 cm) single wheel stem caster set, full swivel - two with brake
- Ergonomic side push handles can be mounted on either side
- Full extension drawers slides for easy access
- Body has pre-threaded accessory mounting holes
- Durable scratch resistant heat sealed powder coat finish
- Fully customizable drawer configurations



Drawer Configuration:

- 18.5" W x 17" D (47 cm x 43.2 cm)
- Three 3" drawers (7.6 cm), one 6" drawer (15.2 cm) and one 9" drawer (22.7 cm)

Warranty: 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts

All dimensions are available in the "MA-Dims-24wide-3DScaster.pdf" document. Shipping weights can be found at the online product listing for this model, or feel free to request them from your Harloff Team.



View similar products online at

<https://www.harloff.com/product-category/emergency-crash-carts/>

Toll Free: (800) 433-4064

medicalsales@harloff.com

www.harloff.com • Ph: (719) 637-0300