



# M-series

## M3DS2424B05

### M-SERIES SHORT EMERGENCY CART with Five Drawers and a Breakaway Lock

#### Features:

- Sturdy, uni-body steel construction
- Choice of lock options
- Heat sealed powder coat steel breakaway lock bar
- Pack of 100 breakaway seals
- Durable scratch resistant heat sealed powder coat finish
- Full extension drawers slides for easy access
- Body has pre-threaded accessory mounting holes
- Removable plastic top
- Ergonomic side push handles can be mounted on either side
- 24" of vertical drawer space
- Drawer Configuration: Three 3" drawers (7.6 cm), one 6" drawer (15.2 cm) and one 9" drawer (22.8 cm)
- Cart is delivered assembled
- Fully customizable drawer configurations
- 3" (7.6 cm) single wheel stem caster set, full swivel - two with brake
- Latex free
- 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts

All dimensions are available in the "MA-Dims-24wide-3DScaster.pdf" document. Shipping weights can be found at the online product listing for this model, or feel free to request them from your Harloff Team.

**Color shown is Red. More than 20 additional colors are available.**



#### Short Height, Medium Width, Five Drawer Emergency Cart with Breakaway Lock

Designed and built to meet the needs of today's healthcare market. Ergonomic side handles, comfort drawer pulls, expanded lock options and a standard plastic top make things easier for your staff and help support your infection control initiatives. M-Series carts offer you choices.



**Engineered &  
Manufactured  
in the U.S.**

**HARLOFF**  
INNOVATING STORAGE FOR CARE

View similar products online at  
<https://www.harloff.com/product-category/emergency-crash-carts/>

Toll Free: (800) 433-4064 • [medicalsales@harloff.com](mailto:medicalsales@harloff.com)  
[www.harloff.com](http://www.harloff.com) • Ph: (719) 637-0300