



M-series

M3DS2430E06+MD24-ANS

M-SERIES TALL EMERGENCY CART with MD24-EMG Emergency Package

Tall, Medium Width, Six Drawer Anesthesia Cart with Electronic Pushbutton Lock

Designed and built to meet the needs of today's healthcare market. Ergonomic side handles, comfort drawer pulls, expanded lock options and a standard antimicrobial top make things easier for your staff and help support your infection control initiatives. M-Series carts offer you choices.



Engineered &
Manufactured
in the U.S.

MD30-ANS Accessories Include:



3-Tier Back Rail
(MD24-RBR3)



(2) 12" Side Rails
(SIDERAIL12)



5-Compartment
Tilt Bin w/Rail Clips
(TILTBIN5RC)



Color shown is Navy.
More than 20 additional
colors are available.

Features:

- Sturdy, uni-body steel construction
- Basic electronic pushbutton lock with key lock override
- Factory set user and supervisor code
- Cabinet may be programmed in the field
- 30" of vertical drawer space
- Cart is delivered assembled
- Antimicrobial plastic top, removable
- Latex free
- 3" (7.6 cm) single wheel stem caster set, full swivel - two with brake
- Ergonomic side push handles can be mounted on either side
- Full extension drawers slides for easy access
- Body has pre-threaded accessory mounting holes
- Durable scratch resistant heat sealed powder coat finish
- Fully customizable drawer configurations



View similar products online at

<https://www.harloff.com/product-category/anesthesia-carts/>

Toll Free: (800) 433-4064

medicalsales@harloff.com

www.harloff.com • Ph: (719) 637-0300

Drawer Configuration:

- 18.5" W x 17" D (47 cm x 43.2 cm)
- Three 3" drawers (7.6 cm), two 6" drawers (15.2 cm) and one 9" drawer (22.7 cm)

Assembled Cart Dimensions:

- 40.75" H x 30" W x 22" D (103.5 cm H x 76.2 cm W x 55.9 cm D) cart only,
- 64.8" H x 32.7" W x 22" D (164.6 cm H x 83.1 cm W x 55.9 cm D) w/accessories
- Cart net weight: 111 pounds (50.3 kg)
- Approx shipping weight: 150 pounds (68 kg)
- Approx shipping dims: 50" H x 33" W x 26" D (127 cm x 99 cm x 66 cm)

Warranty: 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts