



M-series *Drawer Alternative*

EQUIPMENT HOLDER ACCESSORY

for 18", 24" and 30" Wide M-Series Carts

Replace 9", 12 or 18" of vertical drawer space at the top or bottom of your M-Series cart with an open equipment holder, perfect for storing electronics or large medical equipment.

- 17" deep open equipment holding area replaces standard drawer areas in three height options
- 2" pass-through hole in the rear for cords and cables

Dimensions	MDS30 body	MDS24 body	MDS18 body
9" vertical drawer space	MD30-DRW9-EH 8" H x 24.625" W	MD24-DRW9-EH 8" H x 18.625" W	MD18-DRW9-EH 8" H x 12.625" W
12" vertical drawer space	MD30-DRW12-EH 11.125" H x 24.625" W	MD24-DRW12-EH 11.125" H x 18.625" W	MD18-DRW12-EH 11.125" H x 12.625" W
18" vertical drawer space	MD30-DRW18-EH 17.375" H x 24.625" W	MD24-DRW18-EH 17.375" H x 18.625" W	MD18-DRW18-EH 17.375" H x 12.625" W

M-Series Features:

- Sturdy, uni-body steel construction
- Choice of lock options:
www.harloff.com/locking-options/
- Up to 30" of fully customizable vertical space for both drawer and equipment holder configurations
- Full extension drawer slides for easy access
- Cart is delivered assembled
- Antimicrobial plastic top, removable
- Industry only fully curved protective bumper
- Durable scratch resistant hard baked powder coat finish
- Body has pre-threaded accessory mounting holes
- Industry best 5" quiet rolling Tente® casters - 1 steering caster and 2 braking casters
- Latex free
- 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments.
- 1-Year Limited Warranty on any bought-out parts

Colors:

Shown in Red and Forest Green, but available in more than twenty color options which you can view online at www.harloff.com/finish-options/.



View similar products online at <https://www.harloff.com/m-series-carts/>

Toll Free: (800) 433-4064 • medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300