



# M-series

## MDS1821E14

### M-SERIES X-SHORT MINI WIDTH CART

Four Drawers with an Electronic Pushbutton Lock

#### M-Series Features:

- Sturdy, uni-body steel construction
- Basic electronic pushbutton lock
- Key lock override
- Factory set user and supervisor code
- Cabinet may be programmed in the field
- Durable scratch resistant heat sealed powder coat finish
- Full extension drawers slides for easy access
- Body has pre-threaded accessory mounting holes
- Antimicrobial plastic top, removable
- 21" of vertical drawer space
- Cart is delivered assembled
- Industry only fully curved protective bumper
- Fully customizable drawer configurations
- Industry best 5" quiet rolling Tente® casters
  - 1 steering caster and 2 braking casters
- Latex free

#### Warranty:

- 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments.
- 1-Year Limited Warranty on any bought-out parts

#### Colors:

Shown with a Beige body and Light Gray drawers, but available in more than twenty color options which you can view online at [www.harloff.com/finish-options/](http://www.harloff.com/finish-options/).



Engineered & Manufactured in the U.S.



#### Overall Cart Dimensions:

- Internal Drawer Dimensions: 12.5" W x 17" D (31.8 cm W x 43.2 cm D)
- Work Surface Height: 32.875" (83.5 cm)
- Cart dims: 34" H x 23.875" W x 22" D (86.4 cm H x 60.6cm W x 55.9cm D)
- Cart net weight: 81 lbs (36.7 kg)
- Approx shipping weight: 106 lbs (48.1 kg)
- Approx shipping dims: 41" H x 29" W x 24" D (104 cm H x 73.7 cm W x 61 cm D)

#### Drawer Configuration:

- Three 3" drawers (7.6 cm), one 6" drawer (15.2 cm) and one 9" drawer (22.8 cm)

Designed and built to meet the needs of today's healthcare market. Comfort drawer pulls, expanded lock options and a standard antimicrobial top make things easier for your staff and help support your infection control initiatives. M-Series carts offer you choices.

**HARLOFF**  
INNOVATING STORAGE FOR CARE

View similar products online at  
<https://www.harloff.com/m-series-carts/>

Toll Free: (800) 433-4064 • [medicalsales@harloff.com](mailto:medicalsales@harloff.com)  
[www.harloff.com](http://www.harloff.com) • Ph: (719) 637-0300