



**M-series**

**MDS1824B05+MD18-EMG**

## M-SERIES SHORT EMERGENCY CART with MD18-EMG Emergency Package

### Short Mini Width, Five Drawer Emergency Cart w/Breakaway Lock

Designed and built to meet the needs of today's healthcare market. Ergonomic side handles, comfort drawer pulls, expanded lock options and a standard plastic top make things easier for your staff and help support your infection control initiatives. M-Series carts offer you choices.

### MD18-EMG Accessories Include:



Oxygen tank holder (O2HLDR)



Cardiac board and brackets (MD-CARDBRD)



Aluminum 16" Push Side Rail (SIDERAIL16)



IV pole (MD-IV2)



Engineered & Manufactured  
in the U.S.

### Features:

- Sturdy, uni-body steel construction
- Heat sealed powder coat steel breakaway lock bar
- Pack of 100 breakaway seals
- Durable scratch resistant heat sealed powder coat finish
- Full extension drawers slides for easy access
- Body has pre-threaded accessory mounting holes
- Removable plastic top
- 24" of vertical drawer space
- Drawer Configuration: Three 3" drawers (7.6 cm), one 6" drawer (15.2 cm) and one 9" drawer (22.8 cm)
- Cart is delivered assembled
- Industry only fully curved protective bumper
- Fully customizable drawer configurations
- Industry best 5" quiet rolling casters - one steering caster and two braking casters
- 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts

Color shown is Red.  
More than 20 additional  
colors are available.



**HARLOFF**  
INNOVATING STORAGE FOR CARE

View similar products online at

<https://www.harloff.com/product-category/emergency-crash-carts/>

**Toll Free: (800) 433-4064 • Ph: (719) 637-0300**  
**medicalsales@harloff.com • www.harloff.com**

All dimensions are available in the "MA-Dims-18wide-5caster.pdf" document. Shipping weights can be found at the online product listing for this model, or feel free to request them from your Harloff Team.