



M-series

MDS1824B05+MD18-EMG

M-SERIES SHORT EMERGENCY CART with MD18-EMG Emergency Package

Short Mini Width, Five Drawer Emergency Cart w/Breakaway Lock

Designed and built to meet the needs of today's healthcare market. Ergonomic side handles, comfort drawer pulls, expanded lock options and a standard antimicrobial top make things easier for your staff and help support your infection control initiatives. M-Series carts offer you choices.

MD18-EMG Accessories Include:



Oxygen tank holder (O2HLDR)



Cardiac board and brackets (MD-CARDBRD)



Aluminum 16" Push Side Rail (SIDERAIL16)



IV pole (MD-IV2)



Engineered & Manufactured in the U.S.

Features:

- Sturdy, uni-body steel construction
- Heat sealed powder coat steel break-away lock bar
- Pack of 100 break-away seals
- Durable scratch resistant heat sealed powder coat finish
- Full extension drawers slides for easy access
- Body has pre-threaded accessory mounting holes
- Antimicrobial plastic top, removable
- 24" of vertical drawer space
- Cart is delivered assembled
- Industry only fully curved protective bumper
- Fully customizable drawer configurations
- Industry best 5" quiet rolling Tente® casters - 1 steering caster and 2 braking casters

Warranty:

- 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts

Color shown is Red. More than 20 additional colors are available.



View similar products online at

<https://www.harloff.com/product-category/emergency-crash-carts/>

Toll Free: (800) 433-4064 • Ph: (719) 637-0300
medicalsales@harloff.com • www.harloff.com

Drawer Configuration:

- Three 3" drawers (7.6 cm), one 6" drawer (15.2 cm) and one 9" drawer (22.8 cm)
- 12.5" W x 17" D (31.8 cm W x 43.2 cm D)

Assembled Cart Dimensions (without accessories):

- Cart dims: 36.38" H x 23.13" W x 22" D (92.4 cm H x 58.7 cm W x 55.9 cm D)
- Cart net weight: 86 lbs (39 kg)
- Approx shipping weight: 120 lbs (54.4 kg)
- Approx shipping dims: 50" H x 24" W x 24" D (127 cm H x 61 cm W x 61 cm D)