



M-series

MDS1824EKC05

M-SERIES SHORT MINI WIDTH CART with Five Drawers and an Electronic Lock w/Keypad

M-Series Features:

- Sturdy, uni-body steel construction
- Electronic locking cart with keypad access and Key Lock override
- Drawers automatically re-lock after access (pre-set or programmed time)
- Access rights may be assigned for up to 250 users per cart
- May add, edit, or delete user and user groups
- Battery operated with six D-cell Batteries; 1 year minimum life at 25-50 uses per day
- Lock has an LED low battery indicator
- Durable scratch resistant heat sealed powder coat finish
- Full extension drawers slides for easy access
- Body has pre-threaded accessory mounting holes
- Antimicrobial plastic top, removable
- 24" of vertical drawer space
- Cart is delivered assembled
- Industry only fully curved protective bumper
- Fully customizable drawer configurations
- Industry best 5" quiet rolling Tente® casters - 1 steering caster and 2 braking casters
- Latex free

Warranty:

- 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments.
- 1-Year Limited Warranty on any bought-out parts



Engineered & Manufactured in the U.S.



Color shown is White with Navy drawers. More than 20 additional colors are available.

Designed and built to meet the needs of today's healthcare market. Comfort drawer pulls, expanded lock options and a standard antimicrobial top make things easier for your staff and help support your infection control initiatives. M-Series carts offer you choices.

Overall Cart Dimensions:

- Internal Drawer Dimensions: 12.5" W x 17" D (31.8 cm W x 43.2 cm D)
- Work Surface Height: 36.125" (91.8 cm)
- Cart dims: 37.25" H x 23.875" W x 22" D (94.6 cm H x 60.6cm W x 55.9cm D)
- Cart net weight: 86 lbs (39 kg)
- Approx shipping weight: 120 lbs (54.4 kg)
- Approx shipping dims: 50" H x 24" W x 24" D (127 cm H x 61 cm W x 61 cm D)

Drawer Configuration:

- Three 3" drawers (7.6 cm), one 6" drawer (15.2 cm) and one 9" drawer (22.8 cm)

HARLOFF

INNOVATING STORAGE FOR CARE

View similar products online at
<https://www.harloff.com/m-series-carts/>

Toll Free: (800) 433-4064 • medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300