



**M-series**

**MDS1830E06+MD18-ANS**

## M-SERIES TALL ANESTHESIA CART

### Mini Six Drawer Cart w/MD18-ANS Accessory Package

### Tall Height, Mini Width, Six Drawer Anesthesia Cart with an Electronic Pushbutton Lock

Designed and built to meet the needs of today's healthcare market. Comfort drawer pulls, expanded lock options and a standard plastic top make things easier for your staff and help support your infection control initiatives. M-Series carts offer you choices.

### MD18-ANS Accessories Include:



3-Tier Back Rail  
(MD18-RBR3)



(3) 12" Side Rails  
(SIDERAIL12)



5-Compartment  
Tilt Bin w/Rail Clips  
(TILTBIN5RC)



Engineered &  
Manufactured  
in the U.S.

### Features:

- Sturdy, uni-body steel construction
- Basic electronic pushbutton lock with key lock override
- Factory set user and supervisor code
- Cabinet may be programmed in the field
- Durable scratch resistant heat sealed powder coat finish
- Full extension drawers slides for easy access
- Body has pre-threaded accessory mounting holes
- Removable plastic top
- Cart is delivered assembled
- Fully customizable drawer configs with 30" of vertical drawer space
- Industry best 5" quiet rolling casters - one steering caster and two braking casters

### Warranty:

- 12-Year Manufacturer's Warranty on shipments to US and Canada. 1-Year Limited Warranty on any bought-out parts.

Color shown  
with a Light  
Gray body and  
Hammertone  
Blue drawers.  
More than 20  
additional colors  
are available.



**HARLOFF**  
INNOVATING STORAGE FOR CARE

View similar products online at

<https://www.harloff.com/product-category/anesthesia-carts/>

**Toll Free: (800) 433-4064 • Ph: (719) 637-0300**  
**medicalsales@harloff.com • www.harloff.com**

#### Drawer Configuration:

- Three 3" drawers (7.6 cm), two 6" drawers (15.2 cm) & one 9" drawer (22.8 cm)

All dimensions are available in the "MA-Dims-18wide-5caster.pdf" document. Shipping weights can be found at the online product listing for this model, or feel free to request them from your Harloff Team.