



Wood Style

MDS1830EC4-AN

M-SERIES TALL WOOD STYLE CART w/Four Drawers and an Electronic Pushbutton Lock

Why Wood Style Drawer Fronts?

Medical carts which appear more like furniture tend to provide better psychological benefits for your patients, as they feel more at home and comfortable than they would in a sterile setting.

Features:

- **Wood style drawer fronts giving the cart a furniture feel are available in multiple finish options**
- Sturdy, uni-body steel construction
- Basic electronic pushbutton lock
- Key lock override
- Factory set user and supervisor code
- Cabinet may be programmed in the field
- Durable scratch resistant heat sealed powder coat finish on body
- Full extension drawers slides for easy access
- Pre-threaded mounting holes for accessories found at harloff.com/product-category/accessories/
- Antimicrobial plastic top, removable
- 30" of vertical drawer space
- Cart is delivered assembled
- Industry only fully curved protective bumper
- Fully customizable drawer configurations
- Industry best 5" quiet rolling Tente® casters
 - 1 steering caster and 2 braking casters

Two standard wood vinyl colors available plus additional special order vinyl color options with a minimum cart purchase and extended lead time.



Engineered & Manufactured in the U.S.

Colors:

Body shown in Cream, but available in more than twenty color options which you can view online at www.harloff.com/finish-options/. Drawers shown in Wood Vinyl American Natural color.

Overall Cart Dimensions:

- Internal Drawer Dimensions: 12.5" W x 17" D (31.8 cm W x 43.2 cm D)
- Work Surface Height: 41.5" (105.4 cm)
- Cart dims: 42.625" H x 23.875" W x 22" D (108.3 cm H x 60.6cm W x 55.9cm D)
- Cart net weight: 95 lbs (43.1 kg)
- Approx shipping weight: 120 lbs (54.4 kg)
- Approx shipping dims: 50" H x 29" W x 24" D (127 cm H x 73.7 cm W x 61 cm D)

Drawer Configuration:

- One 3" drawer (7.6 cm) and three 9" drawers (22.8 cm)

Warranty: 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts



Hospitals, Surgery Centers and Clinics can bring the warm appearance of wood to their treatment and other procedure carts with Harloff's line of M-Series carts with vinyl-coated aluminum drawers. Designed and built to combine the aesthetics of wood with the long-term functionality of a Harloff metal cart.



View similar products online at

<https://www.harloff.com/product-category/m-series-medical-carts/>

Toll Free: (800) 433-4064 • medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300