



M-series

MDS1830K06+MD18-PHB

M-SERIES TALL PHLEBOTOMY CART

Mini Six Drawer Cart w/MD18-PHB Accessory Package

Tall Height, Mini Width, Six Drawer Phlebotomy Cart with Standard Key Lock

Designed and built to meet the needs of today's healthcare market. Phlebotomy carts keep your supplies nearby, are mobile, provide a dedicated work surface for procedures and allow quick disposal of used needles. M-Series carts offer you choices.

MD18-PHB Accessories Include:



6-Compartment
Tilt Bin w/Rail Clips
(TILTBIN6RC)



Glove Box Holder
Single Direct Mount
GLOVEBOX1DM



Aluminum Waste
Container with Lid
AL680436

Engineered &
Manufactured
in the U.S.



Push Handle
(MINI-PBHANDKIT)



3-Tier Back Rail
(MD18-RBR3)

Color shown
with a Cream
body and
White drawers.
More than 20
additional colors
are available.



- Sturdy, uni-body steel construction
- Key lock with two keys
- Durable scratch resistant heat sealed powder coat finish
- Full extension drawer slides for easy access
- Body has pre-threaded accessory mounting holes
- Removable plastic top
- Cart is delivered assembled
- Industry best 5" quiet rolling casters - one steering caster and two braking casters
- Fully customizable drawer configs with 30" of vertical drawer space

HARLOFF
INNOVATING STORAGE FOR CARE

View similar products online at

<https://www.harloff.com/product-category/phlebotomy-carts-iv-start-carts/>

Toll Free: (800) 433-4064 • Ph: (719) 637-0300
medicalsales@harloff.com • www.harloff.com

Drawer Configuration:

- Three 3" drawers (7.6 cm), two 6" drawers (15.2 cm) & one 9" drawer (22.8 cm)

Warranty: 12-Year Manufacturer's Warranty on shipments to US and Canada.
1-Year Limited Warranty on any bought-out parts.

All dimensions are available in the "MA-Dims-18wide-5caster.pdf" document. Shipping weights can be found at the online product listing for this model, or feel free to request them from your Harloff Team.