



MALIGNANT HYPERTHERMIA

MDS2430B02-MH1

MALIGNANT HYPERTHERMIA CART

with 1.0 Cubic Feet Medical Grade Refrigerator

Features:

- Low carbon powder-coated steel construction w/built-in Accucold® ARS1PV refrigerator
- Cart holds critical MH emergency supplies in one place and at the proper temperature in close proximity to areas where most needed
- High quality Accucold® refrigerator ensures temperature control for insulin, cold packs and cold saline solution for IV infusion and/or gastric, peritoneal or rectal irrigation, as appropriate.
- Secured drawers for storage of injectables, syringes, catheters, NG tubes, filters, and other necessary supplies
- Breakaway locks on drawers and refrigerator (refrigerator can be keylocked) with 100 breakaway seals
- Full extension drawer slides for easy access
- Industry only fully curved protective bumper
- Industry best 5" quiet rolling casters - one steering caster and two braking casters
- Cart ships ready to use with wheels installed
- Removable plastic top

M-series

Colors:

Shown in Navy, but available in more than twenty color options. View them online: www.harloff.com/finish-options/



**Malignant Hyperthermia
Cart, Two Drawer, Breakaway
Lock with No Accessories**

**Engineered & Manufactured
in the U.S.**

Warranty:

- 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments.
- 1-Year Limited Warranty on any bought-out parts

Drawer Configuration:

- One 3" (7.6 cm) drawer, one 6" (15.2 cm) drawer & (1) Compartment with pre-installed refrigerator

Optional Accessories:



Large wire basket



Side storage
cabinets



Two large chart/
document holders



Waste container &
mounting bracket

All dimensions are available in the "MA-Dims-24wide-5caster.pdf" document. Shipping weights can be found at the online product listing for this model, or feel free to request them from your Harloff Team.

HARLOFF
INNOVATING STORAGE FOR CARE

View similar products online at

<https://www.harloff.com/product-category/malignant-hyperthermia-carts/>

**Toll Free: (800) 433-4064 • medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300**