



# M-series

## MDS2430E06

### M-SERIES MEDIUM WIDTH CART w/Six Drawers and Electronic Pushbutton Lock

#### Features:

- Sturdy, uni-body steel construction
- Basic electronic pushbutton lock
- Key lock override
- Factory set user and supervisor code
- Cabinet may be programmed in the field
- Durable scratch resistant heat sealed powder coat finish
- Full extension drawers slides for easy access
- Body has pre-threaded accessory mounting holes
- Antimicrobial plastic top, removable
- Ergonomic side push handles can be mounted on either side
- 30" of vertical drawer space
- Cart is delivered assembled
- Industry only fully curved protective bumper
- Fully customizable drawer configurations
- Industry best 5" quiet rolling Tente® casters - 1 steering caster and 2 braking casters
- Latex free



Engineered & Manufactured in the U.S.



Color shown is Light Gray with Hammertone Blue drawers. More than 20 additional colors are available.

#### Tall, Medium Width, Six Drawer Cart with Electronic Pushbutton Lock

Designed and built to meet the needs of today's healthcare market. Ergonomic side handles, comfort drawer pulls, expanded lock options and a standard antimicrobial top make things easier for your staff and help support your infection control initiatives. M-Series carts offer you choices.

#### Drawer Configuration:

- 18.5" W x 17" D (47 cm x 43.2 cm)
- Three 3" drawers (7.6 cm), Two 6" drawers (15.2 cm) and one 9" drawer (22.9 cm)

#### Assembled Cart Dimensions:

- 43.75" H x 30" W x 22" D (111.1 cm H x 76.2 cm W x 55.9 cm D)
- Cart net weight: 119 pounds (54 kg)
- Approx shipping weight: 158 pounds (71.7 kg)
- Approx shipping dims: 56" H x 33" W x 26" D (142.2 cm x 99 cm x 66 cm)

**Warranty:** 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts

# HARLOFF

INNOVATING STORAGE FOR CARE

View similar products online at

<https://www.harloff.com/product-category/m-series-medical-carts/>

Toll Free: (800) 433-4064 • [medicalsales@harloff.com](mailto:medicalsales@harloff.com)  
[www.harloff.com](http://www.harloff.com) • Ph: (719) 637-0300