



# M-series

## MDS2430K06+MD24-ANS

### M-SERIES TALL ANESTHESIA CART Six Drawer Cart with MD24-ANS Accessory Package

#### Tall Height, Medium Width, Six Drawer Anesthesia Cart with a Standard Key Lock

Designed and built to meet the needs of today's healthcare market. Comfort drawer pulls, expanded lock options and a standard antimicrobial top make things easier for your staff and help support your infection control initiatives. M-Series carts offer you choices.

#### MD24-ANS Accessories Include:



3-Tier Back Rail (MD24-RBR3)



(3) 12" Side Rails (SIDERAIL12)



5-Compartment Tilt Bin w/Rail Clips (TILTBIN5RC)



Engineered & Manufactured in the U.S.



Color shown with a **Hammeritone Blue** body and **Hammeritone Gray** drawers. More than 20 additional colors are available.

#### Features:

- Sturdy, uni-body steel construction with key lock and two keys
- Durable scratch resistant heat sealed powder coat finish
- Full extension drawers slides for easy access
- Antimicrobial plastic top, removable
- Ergonomic side push handles can be mounted on either side
- Body has pre-threaded accessory mounting holes
- 30" of vertical drawer space
- Cart is delivered assembled
- Industry only fully curved protective bumper
- Fully customizable drawer configurations
- Industry best 5" quiet rolling Tente® casters
  - 1 steering caster and 2 braking casters

#### Warranty:

- 12-Year Manufacturer's Warranty on shipments to US and Canada. 1-Year Limited Warranty on any bought-out parts.



View similar products online at

<https://www.harloff.com/product-category/anesthesia-carts/>

**Toll Free: (800) 433-4064 • Ph: (719) 637-0300**  
**medicalsales@harloff.com • www.harloff.com**

#### Overall Cart Dimensions:

- Work Surface Height: 41.5" (105.4 cm)
- Internal Drawer Dimensions: 12.5" W x 17" D (31.8 cm W x 43.2 cm D)
- Cart dims: 42.625" H x 23.875" W x 22" D (108.3 cm H x 60.6cm W x 55.9cm D)
- Cart net weight: 95 lbs (43.1 kg)
- Approx shipping weight: 120 lbs (54.4 kg)
- Approx shipping dims: 50" H x 29" W x 24" D (127 cm H x 73.7 cm W x 61 cm D)

#### Drawer Configuration:

- Three 3" drawers (7.6 cm), two 6" drawers (15.2 cm) & one 9" drawer (22.8 cm)