



M-series

MDS2430KC5 with Equipment Holder

M-SERIES MEDIUM WIDTH CART with 3 Drawers, 18" Equipment Holder and Key Lock

M-Series Features:

- Sturdy, uni-body steel construction
- Choice of lock options:
www.harloff.com/locking-options/
- Up to 30" of fully customizable vertical space is available for both drawer and equipment holder configurations
- 17" deep open equipment holding area can replace standard drawer areas in three possible height options
- 2" pass-through hole in the rear for cords and cables
- Full extension drawer slides for easy access
- Cart is delivered assembled
- Antimicrobial plastic top, removable
- Industry only fully curved protective bumper
- Durable scratch resistant hard baked powder coat finish
- Body has pre-threaded accessory mounting holes
- Industry best 5" quiet rolling Tente® casters - 1 steering caster and 2 braking casters
- Latex free

Warranty:

- 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments. 1-Year Limited Warranty on any bought-out parts

Colors:

Shown in Forest Green, but available in more than twenty color options which you can view online at www.harloff.com/finish-options/.

Overall Cart Dimensions:

- Internal Drawer Dimensions: 18.5" W x 17" D (47 x 43.2 cm)
- Work Surface Height: 42.625" (108.3 cm)
- Cart dims: 43.75" H x 30" W x 22" D (111.1 cm H x 76.2cm W x 55.9cm D)
- Cart net weight: 119 pounds (54 kg)
- Approx shipping weight: 158 pounds (71.7 kg)
- Approx shipping dims: 56" H x 33" W x 26" D (142.2 cm x 99 cm x 66 cm)

Drawer Configuration:

- Two 3" drawers (7.6 cm), one 6" drawer (15.2 cm) and one MD24-DRW18-EH equipment holder with an 17.4" H x 18.63" W x 17" D (44.1 cm H x 47.3cm W x 43.2cm D) opening



Designed and built to meet the needs of today's healthcare market. Replace 9", 12 or 18" of vertical drawer space at the top or bottom of your M-Series cart with an open equipment holder, perfect for storing electronics or large medical equipment. M-Series carts offer you choices.

HARLOFF

INNOVATING STORAGE FOR CARE

View similar products online at
<https://www.harloff.com/m-series-carts/>

Toll Free: (800) 433-4064 • medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300