



M-SERIES X-SHORT EMERGENCY CART with Five Drawers and a Breakaway Lock

•

HARLOFF

Features:

- · Sturdy, uni-body steel construction
- Choice of lock options
- Heat sealed powder coat steel break-away lock bar
- Pack of 100 break-away seals
- · Durable scratch resistant heat sealed powder coat finish
- · Full extension drawers slides for easy access
- Body has pre-threaded accessory mounting holes
- · Antimicrobial plastic top, removable
- · Ergonomic side push handles can be mounted on either side
- 21" of vertical drawer space
- Cart is delivered assembled
- · Industry only fully curved protective bumper
- Fully customizable drawer configurations
- Industry best 5" quiet rolling Tente® casters - 1 steering caster and 2 braking casters
- Latex free

Color shown is Red. More than 20 additional colors are available.

Engineered & Manufactured in the U.S.

X-Short Height, Standard Width, **Five Drawer Emergency Cart** with Breakaway Lock

Designed and built to meet the needs of today's healthcare market. Ergonomic side handles, comfort drawer pulls, expanded lock options and a standard antimicrobial top make things easier for your staff and help support your infection control initiatives. M-Series carts offer you choices.

Drawer Configuration:

- Four 3" drawers (7.6 cm) and one 9" drawer (22.9 cm)
- 24.5" W x 17" D (62.2 x 43.2 cm)

Assembled Cart Dimensions:

- 34" H x 36.75" W x 22" D (86.4 cm H x 93.3 cm W x 55.9 cm D)
- Cart weight: 129 lbs (58.5 kg)
- Approximate shipping weight: 168 lbs (76.2 kg)
- Approximate Shipping Dims: 56" H x 37" W x 26" D (142.2 cm x 94 cm x 66 cm)

Warranty:

12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments. 1-Year Limited Warranty on any bought-out parts



View similar products online at

https://www.harloff.com/acute-care-carts/emergency-crash-carts/

Toll Free: (800) 433-4064 • medicalsales@harloff.com www.harloff.com • Ph: (719) 637-0300