



M-series MDS3024B05+MD30-EMG1

M-SERIES SHORT EMERGENCY CART with MD30-EMG1 Emergency Package

Short, Standard Width, 5-Drawer Emergency Cart w/Breakaway Lock

Designed and built to meet the needs of today's healthcare market. Ergonomic side handles, comfort drawer pulls, expanded lock options and a standard antimicrobial top make things easier for your staff and help support your infection control initiatives. M-Series carts offer you choices.

MD30-EMG1 Accessories Include:



Drop shelf (MD-DRPSHLF)



Oxygen tank holder (O2HLDR)



IV pole (MD-IV2)

Color shown is Red. More than 20 additional colors are available.



Cardiac board and brackets (MD-CARDBRD)



(2 Sets) Standard 3" Drawer Dividers (MD30-DIV3-S)



(2) Utility hooks (UHOOKS2)



Features:

- Sturdy, uni-body steel construction
- Choice of lock options
- Pack of 100 break-away seals
- 24" of vertical drawer space
- Full extension drawers slides for easy access
- Cart is delivered assembled
- Antimicrobial plastic top, removable
- Latex free
- Fully customizable drawer configurations
- Industry best 5" quiet rolling Tente® casters - 1 steering caster and 2 braking casters
- Heat sealed powder coat steel break-away lock bar
- Ergonomic side push handles can be mounted on either side
- Body has pre-threaded accessory mounting holes
- Industry only fully curved protective bumper
- Durable scratch resistant heat sealed powder coat finish



Engineered & Manufactured in the U.S.

Drawer Configuration:

- Three 3" drawers (7.6 cm), one 6" drawer (15.2 cm) and one 9" drawer (22.9 cm)
- 24.5" W x 17" D (62.2 x 43.2 cm)

Assembled Cart Dimensions:

- 37.25" H x 36.75" W x 22" D (94.6 cm H x 93.3 cm W x 55.9 cm D) cart only,
- 65.31" H x 39.44" W x 22" D (165.9 cm H x 100.2 cm W x 55.9 cm D) w/accessories
- Cart weight: 137 lbs (62.1 kg)
- Approximate shipping weight: 161 lbs (73 kg)
- Approximate Shipping Dims: 56" H x 37" W x 26" D (142.2 cm x 94 cm x 66 cm)

Warranty: 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts



View similar products online at <https://www.harloff.com/product-category/emergency-crash-carts/>

Toll Free: (800) 433-4064 • Ph: (719) 637-0300
medicalsales@harloff.com • www.harloff.com