



**M-series**

**MDS3024E14+MD30-ANS**

## M-SERIES SHORT ANESTHESIA CART

### Four Drawer Cart with MD30-ANS Accessory Package

### Short Height, Standard Width, Four Drawer Cart with an Electronic Pushbutton Lock

Designed and built to meet the needs of today's healthcare market. Ergonomic side handles, comfort drawer pulls, expanded lock options and a standard antimicrobial top make things easier for your staff and help support your infection control initiatives.

### MD30-ANS Accessories Include:



3-Tier Back Rail  
(MD30-RBR3)



5-Compartment  
Tilt Bin w/Rail Clips  
(TILTBIN5RC)



(3) 12" Side Rails  
(SIDERAIL12)



Engineered &  
Manufactured  
in the U.S.

### Features:

- Sturdy, uni-body steel construction
- Basic electronic pushbutton lock with key lock override
- Factory set user and supervisor code
- Cabinet may be programmed in the field
- Durable scratch resistant heat sealed powder coat finish
- Full extension drawers slides for easy access
- Body has pre-threaded accessory mounting holes
- Antimicrobial plastic top, removable
- Ergonomic side push handles mounted on either side
- 24" of vertical drawer space
- Cart is delivered assembled
- Industry only fully curved protective bumper
- Fully customizable drawer configurations
- Industry best 5" quiet rolling Tente® casters - 1 steering caster and 2 braking casters
- Latex free

Color shown in Light Blue.  
More than 20 additional  
colors are available.



**HARLOFF**  
INNOVATING STORAGE FOR CARE

View similar products online at

<https://www.harloff.com/product-category/anesthesia-carts/>

**Toll Free: (800) 433-4064 • Ph: (719) 637-0300**  
**medicalsales@harloff.com • www.harloff.com**

#### Overall Cart Dimensions:

- Cart weight: 137 lbs (62.1 kg)
- 37.25" H x 36.75" W x 22" D (94.6 cm H x 93.3 cm W x 55.9 cm D)
- Approximate shipping weight: 176 lbs (79.8 kg)
- Approximate Shipping Dims: 56" H x 37" W x 26" D (142.2 cm x 94 cm x 66 cm)

#### Drawer Configuration:

- Internal Drawer Dimensions: 24.5" W x 17" D (62.2 x 43.2 cm)
- Two 3" drawers (7.6 cm), one 6" drawer (15.2 cm) and one 12" drawer (30.5 cm)

**Warranty:** 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts