



# M-series

## MDS3024K05

### M-SERIES SHORT MEDICAL CART with Five Drawers and a Standard Key Lock

#### Features:

- Sturdy, uni-body steel construction
- Key lock with two keys
- Durable scratch resistant heat sealed powder coat finish
- Full extension drawer slides for easy access
- Body has pre-threaded accessory mounting holes
- Antimicrobial plastic top, removable
- Ergonomic side push handles can be mounted on either side
- 24" of vertical drawer space
- Cart is delivered assembled
- Industry only fully curved protective bumper
- Fully customizable drawer configurations
- Industry best 5" quiet rolling Tente® casters - 1 steering caster and 2 braking casters
- Latex free

Engineered &  
Manufactured  
in the U.S.



Color shown is a Beige body with White drawers. More than 20 additional colors are available.



#### Short Height, Standard Width, Five Drawer Medical Cart with Standard Key Lock

Designed and built to meet the needs of today's healthcare market. Ergonomic side handles, comfort drawer pulls, expanded lock options and a standard antimicrobial top make things easier for your staff and help support your infection control initiatives. M-Series carts offer you choices.

#### Drawer Configuration:

- Three 3" drawers (7.6 cm), one 6" drawer (15.2 cm) and one 9" drawer (22.9 cm)
- 24.5" W x 17" D (62.2 x 43.2 cm)

#### Assembled Cart Dimensions:

- 37.25" H x 36.75" W x 22" D (94.6 cm H x 93.3 cm W x 55.9 cm D)
- Cart weight: 137 lbs (62.1 kg)
- Approximate shipping weight: 176 lbs (79.8 kg)
- Approximate Shipping Dims: 56" H x 37" W x 26" D (142.2 cm x 94 cm x 66 cm)

**Warranty:** 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts.

**HARLOFF**  
INNOVATING STORAGE FOR CARE

View similar products online at

<https://www.harloff.com/product-category/m-series-medical-carts/>

Toll Free: (800) 433-4064 • [medicalsales@harloff.com](mailto:medicalsales@harloff.com)  
[www.harloff.com](http://www.harloff.com) • Ph: (719) 637-0300