



M-series

MDS3024K14+MD30-ANS

M-SERIES SHORT ANESTHESIA CART

Four Drawer Cart with MD30-ANS Accessory Package

Short Height, Standard Width, Four Drawer Anesthesia Cart with a Standard Key Lock

Designed and built to meet the needs of today's healthcare market. Ergonomic side handles, comfort drawer pulls, expanded lock options and a standard antimicrobial top make things easier for your staff and help support your infection control initiatives. M-Series carts offer you choices.

MD30-ANS Accessories Include:



3-Tier Back Rail
(MD30-RBR3)



5-Compartment
Tilt Bin w/Rail Clips
(TILTBIN5RC)



(3) 12" Side Rails
(SIDERAIL12)

Color shown with a White body and Cream drawers. More than 20 additional colors are available.

Features:

- Sturdy, uni-body steel construction, delivered assembled
- Key lock with two keys
- Durable scratch resistant heat sealed powder coat finish
- Full extension drawers slides for easy access
- Body has pre-threaded accessory mounting holes
- Antimicrobial plastic top, removable
- Ergonomic side push handles, mounted on either side
- 24" of vertical drawer space
- Industry only fully curved protective bumper
- Fully customizable drawer configurations
- Industry best 5" quiet rolling Tente® casters
 - 1 steering caster and 2 braking casters

Warranty:

- 12-Year Manufacturer's Warranty on shipments to US and Canada. 1-Year Limited Warranty on any bought-out parts.



**Engineered &
Manufactured
in the U.S.**

HARLOFF
INNOVATING STORAGE FOR CARE

View similar products online at

<https://www.harloff.com/product-category/anesthesia-carts/>

Toll Free: (800) 433-4064 • Ph: (719) 637-0300
medicalsales@harloff.com • www.harloff.com

Overall Cart Dimensions:

- Internal Drawer Dimensions: 24.5" W x 17" D (62.2 x 43.2 cm)
- 37.25" H x 36.75" W x 22" D (94.6 cm H x 93.3 cm W x 55.9 cm D)
- Cart weight: 137 lbs (62.1 kg)
- Approximate shipping weight: 176 lbs (79.8 kg)
- Approximate Shipping Dims: 56" H x 37" W x 26" D (142.2 cm x 94 cm x 66 cm)

Drawer Configuration:

- Two 3" drawers (7.6 cm), one 6" drawer (15.2 cm) and one 12" drawer (30.5 cm)