



M-series MDS3030B05-12FD

M-SERIES CART WITH FLIP-UP DOOR Crash Cart w/Five Drawers and Breakaway Lock

Includes a lower flip-up door that stows horizontally and, when opened, provides a 12" high inside storage compartment

Features:

- Sturdy, uni-body steel construction
- Choice of lock options
- Heat sealed powder coat steel breakaway lock bar
- Pack of 100 breakaway seals
- Durable scratch resistant heat sealed powder coat finish
- Full extension drawers slides for easy access
- Body has pre-threaded accessory mounting holes
- Removable plastic top
- Ergonomic side push handles can be mounted on either side
- Industry only fully curved protective bumper
- Fully customizable drawer configurations
- Industry best 5" quiet rolling casters - one steering caster and two braking casters
- 30" of vertical drawer space
- Drawer Configuration: Four 3" drawers (7.6 cm), one 6" drawer (15.2 cm) and one 12" flip-up drawer (30.4 cm)
- Cart is delivered assembled
- Latex free
- 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts



**Color shown is Red.
More than 20 additional
colors are available.**

All dimensions are available in the "MA-Dims-30wide-5caster.pdf" document. Shipping weights can be found at the online product listing for this model, or feel free to request them from your Harloff Team.

Tall, Standard Width, Five Drawer Emergency Cart with Flip-Up Door and Breakaway Lock

Designed and built to meet the needs of today's healthcare market. Ergonomic side handles, comfort drawer pulls, expanded lock options and a standard plastic top make things easier for your staff and help support your infection control initiatives. M-Series carts offer you choices.



**Engineered &
Manufactured
in the U.S.**



View similar products online at

<https://www.harloff.com/product-category/emergency-crash-carts/>

**Toll Free: (800) 433-4064 • medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300**