



# M-series

# MDS3030B06-EMG1

## M-SERIES TALL EMERGENCY CART with MD30-EMG1 Emergency Package

### Tall, Standard Width, Six Drawer Emergency Cart w/Breakaway Lock

Designed and built to meet the needs of today's healthcare market. Ergonomic side handles, comfort drawer pulls, expanded lock options and a standard antimicrobial top make things easier for your staff and help support your infection control initiatives. M-Series carts offer you choices.

#### MD30-EMG1 Accessories Include:

Color shown is Red. More than 20 additional colors are available.



Drop shelf  
(MD-DRPSHLF)



Oxygen tank  
holder (O2HLDR)



IV pole (MD-IV2)



Cardiac board  
and brackets  
(MD-CARDBRD)



(2 Sets) Standard  
3" Drawer Dividers  
(MD30-DIV3-S)



(2) Utility hooks  
(UHOOKS2)



#### Features:

- Sturdy, uni-body steel construction
- Choice of lock options
- Pack of 100 break-away seals
- 30" of vertical drawer space
- Full extension drawers slides for easy access
- Cart is delivered assembled
- Antimicrobial plastic top, removable
- Latex free
- Fully customizable drawer configurations
- Industry best 5" quiet rolling Tente® casters - 1 steering caster and 2 braking casters
- Heat sealed powder coat steel break-away lock bar
- Ergonomic side push handles can be mounted on either side
- Body has pre-threaded accessory mounting holes
- Industry only fully curved protective bumper
- Durable scratch resistant heat sealed powder coat finish



Engineered & Manufactured in the U.S.

#### Drawer Configuration:

- 24.5" W x 17" D (62.2 x 43.2 cm)
- Three 3" drawers (7.6 cm), two 6" drawers (15.2 cm) and one 9" drawer (22.7 cm)

#### Assembled Cart Dimensions:

- 43.75" H x 36.75" W x 22" D (111.1 cm H x 93.3 cm W x 55.9 cm D) cart only, 71.81" H x 39.44" W x 22" D (182.4 cm H x 100.2 cm W x 55.9 cm D) w/accessories
- Cart weight: 190 lbs (86.2 kg)
- Approximate shipping weight: 229 lbs (103.9 kg)
- Approximate Shipping Dims: 56" H x 37" W x 26" D (142.2 cm x 94 cm x 66 cm)

**Warranty:** 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments. 1-Year Limited Warranty on any bought-out parts

# HARLOFF

INNOVATING STORAGE FOR CARE

View similar products online at  
<https://www.harloff.com/acute-care-carts/emergency-crash-carts/>

Toll Free: (800) 433-4064 • Ph: (719) 637-0300  
[medicalsales@harloff.com](mailto:medicalsales@harloff.com) • [www.harloff.com](http://www.harloff.com)