



# M-series MDS3030B09PED-EMG

## M-SERIES TALL PEDIATRIC CART with Nine Drawers and Emergency Package

Designed and built to meet the needs of today's healthcare market. Ergonomic side handles, comfort drawer pulls, expanded lock options and a standard antimicrobial top make things easier for your staff and help support your infection control initiatives. M-Series carts offer you choices.

### Accessories Include:



Side Mounted  
Drop Shelf  
(MD-DRPSHLF)



Cardiac board  
and brackets  
(MD-CARDBRD)



Oxygen tank  
holder (O2HLDR)



Utility Hooks, Set  
of 2 (UHOOKS2)



IV pole (MD-IV2)

Engineered &  
Manufactured  
in the U.S.



### Features:

- Sturdy, uni-body steel construction
- Cart is designed to be compatible with popular pediatric color coding systems
- Tall cabinet with separate breakaway locks on each drawer
- Industry only fully curved protective bumper
- Fully customizable drawer configurations
- Durable scratch resistant heat sealed powder coat finish
- Full extension drawers slides for easy access
- Ergonomic side push handles can be mounted on either side
- Comfort drawer pulls
- Body has pre-threaded accessory mounting holes
- Industry best 5" quiet rolling Tente® casters - 1 steering caster & 2 braking casters
- Antimicrobial plastic top, removable
- Painted steel break-away lock bar
- Pack of 100 break-away seals
- 30" of vertical drawer space
- Latex free
- Cart is delivered assembled

**Color shown is Beige with  
color coded drawers. More  
than 20 additional body  
colors are available.**



View similar products online at

<https://www.harloff.com/product-category/pediatric-resuscitation-carts/>

**Toll Free: (800) 433-4064**  
**medicalsales@harloff.com**  
**www.harloff.com • Ph: (719) 637-0300**

#### Drawer Configuration:

- 24.5" W x 17" D (62.2 x 43.2 cm)
- Eight 3" (7.6 cm) drawers and one 6" (15.2 cm) drawer

#### Assembled Cart Dimensions:

- 43.75" H x 36.75" W x 22" D (111.1 cm H x 93.3 cm W x 55.9 cm D)
- Writing surface: 42.625" (108.2 cm)
- Cart weight: 155 lbs (70.5 kg)
- Approximate shipping weight: 195 lbs (88.6 kg)
- Approximate shipping dims: 56" H x 39" W x 26" D (142.2 cm x 99 cm x 66 cm)

**Warranty:** 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments. 1-Year Limited Warranty on any bought-out parts