



# M-series

## MDS3030B08

### M-SERIES TALL EMERGENCY CART with Eight Drawers and a Breakaway Lock

#### Features:

- Sturdy, uni-body steel construction
- Choice of lock options
- Heat sealed powder coat steel break-away lock bar
- Pack of 100 break-away seals
- Durable scratch resistant heat sealed powder coat finish
- Full extension drawers slides for easy access
- Body has pre-threaded accessory mounting holes
- Antimicrobial plastic top, removable
- Ergonomic side push handles can be mounted on either side
- 30" of vertical drawer space
- Cart is delivered assembled
- Industry only fully curved protective bumper
- Fully customizable drawer configurations
- Industry best 5" quiet rolling Tente® casters - 1 steering caster and 2 braking casters
- Latex free



Engineered & Manufactured in the U.S.

Color shown is Light Blue. More than 20 additional colors are available.



#### Tall, Standard Width, Eight Drawer Emergency Cart with Breakaway Lock

Designed and built to meet the needs of today's healthcare market. Ergonomic side handles, comfort drawer pulls, expanded lock options and a standard antimicrobial top make things easier for your staff and help support your infection control initiatives. M-Series carts offer you choices.

#### Drawer Configuration:

- 24.5" W x 17" D (62.2 x 43.2 cm)
- Six 3" drawers (7.6 cm) and two 6" drawers (15.2 cm)

#### Assembled Cart Dimensions:

- 43.75" H x 36.75" W x 22" D (111.1 cm H x 93.3 cm W x 55.9 cm D)
- Cart weight: 150 lbs (68 kg)
- Approximate shipping weight: 189 lbs (85.7 kg)
- Approximate Shipping Dims: 56" H x 37" W x 26" D (142.2 cm x 94 cm x 66 cm)

**Warranty:** 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts



View similar products online at <https://www.harloff.com/product-category/emergency-crash-carts/>

Toll Free: (800) 433-4064 • [medicalsales@harloff.com](mailto:medicalsales@harloff.com)  
[www.harloff.com](http://www.harloff.com) • Ph: (719) 637-0300