



M-series MDS3030K02-21DR

M-SERIES STANDARD WIDTH CART with Two Drawers, Storage Area and Key Lock

Includes a 21" double door storage compartment (keyed separately) and adjustable interior shelf with six possible height locations

Features:

- Sturdy, uni-body steel construction
- Key lock with two keys
- Durable scratch resistant heat sealed powder coat finish
- Full extension drawers slides for easy access
- Body has pre-threaded accessory mounting holes
- Antimicrobial plastic top, removable
- 30" of vertical drawer space
- Cart is delivered assembled
- Industry only fully curved protective bumper
- Industry best 5" quiet rolling Tente® casters - 1 steering caster and 2 braking casters
- Fully customizable drawer configurations
- Latex free



Engineered & Manufactured in the U.S.



Color shown is Navy. More than 20 additional colors are available.

Tall, Standard Width, Two Drawer Cart and Storage Compartment with Key Lock

Designed and built to meet the needs of today's healthcare market. With a side handle, comfort drawer pulls, expanded lock options and a standard antimicrobial top to make things easier for your staff and to help support your infection control initiatives. M-Series carts offer you choices.

Drawer Configuration:

- 24.5" W x 17" D (62.2 x 43.2 cm)
- One 3" drawer (7.6 cm), one 6" drawer (15.2 cm) and one 21" double door storage compartment (53.3 cm)

Assembled Cart Dimensions:

- 43.75" H x 36.75" W x 22" D (111.1 cm H x 93.3 cm W x 55.9 cm D)
- Cart weight: 150 lbs (68 kg)
- Approximate shipping weight: 189 lbs (85.7 kg)
- Approximate Shipping Dims: 56" H x 37" W x 26" D (142.2 cm x 94 cm x 66 cm)

Warranty: 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments. 1-Year Limited Warranty on any bought-out parts



View similar products online at

<https://www.harloff.com/product-category/treatment-and-procedure-carts/>

**Toll Free: (800) 433-4064 • medicalsales@harloff.com
www.harloff.com • Ph: (719) 637-0300**